

The spacer cover promotes heat and moisture dissipation, whilst utilising the padded construction to enhance comfort.

A contoured profile aids pressure distribution while maintaining a high degree of posture control and comfort.

Printed instructions and information.

A concealed zip helps prevent fluid contamination, preferred by infection control officers.

High quality reflex foams.

A low profile contoured laminated insert with built in adduction, abduction and general posture control, whilst maintaining security and comfort.


Durable, heavy duty zip, incorporating a bright yellow slider making it easier to see, locate and use.

Rolled edge cover.

Piping to give an extremely smart appearance.

Non slip base for added security and greater client confidence.

Vapour permeable-thermo plastic membrane gives clinical protection against contamination.



APPLICATION

***PAEDIATRIC CUSHION
WITH POSTURE CONTROL**

Sizes: 9"x9" to 14"x15".

Mercury Junior

PC



MMS Medical
51 East Gate Drive,
Little Island,
Cork.
T: 021 461 8000
F: 021 461 8099
E: sales@mmsmedical.ie

www.mmsmedical.ie

MMS Medical
Unit 81,
Western Parkway Business Park,
Lower Ballymount Road,
Dublin 12.
T: 01 450 0582
F: 01 456 5889
E: sales@mmsmedical.ie



Class 1

FM 30965

Mercury Junior

- Description:** A paediatric positional component, aimed at providing posture control. The Mercury Junior offers long term durability, stability, comfort and support. Available in either flat or crescent based format with the option of enhanced pressure distribution with the addition of memory foam.
- Sizes:** 9"x9" to 14"x15". All cushions are 2 ¼" in thickness.
- Cover:** A spacer cover is issued as standard equipment, chosen for its superior heat and moisture dissipation and bringing with it a high degree of comfort, trimmed with piping and vapour permeable side banding, constructed on a non slip base for added security and incorporating a high visibility quality zip and slider, concealed for increased clinical performance. A vapour permeable thermo plastic interliner is included as standard. All relevant technical, directional and care information is clearly printed on the rear of the main cover.
- Risk Category*:** Low/medium
- Weight Guide†:** For maximum comfort, durability and performance the weight guide on a 14" x 15" is 10 stone (64kg)
- Support Pressure:** 1 lb per square inch at 50% compression or deflection
- Linear Load Limit:** 21 stone (133kg) on a typical 17"x17" cushion

* Risk categories are given as a guide and rely on your judgement and experience.

† Weight Guides are a general indication of an ideal user weight to give optimum durability, comfort and performance; it is not a definitive or maximum value. (This product has been tested to 325lb (146kg) without failure).

Instructions for Use and Cleaning

DESIGN

The device has been designed to provide the correct levels of pressure reduction, support and comfort which are needed and according to the design the correct amount of postural assistance.

PLACEMENT

1. Before using the cushion for the first time, please remove the outer packaging if present and check that the chair seat has no objects or projections which may damage the new cushion.
2. Place the device onto the chair ensuring that our label is on the base of the device (ie. do not sit or lean on the label side).

CLEANING

The cushion has removable covers and can be cleaned according to the material type. If possible sponge off any spillage or contamination as soon as possible with hot water and soap. Do not remove the special transparent liner (if fitted) but sponge it down.

DO NOT PLACE THE WHOLE DEVICE DIRECTLY INTO A WASHING MACHINE

DOMESTIC CLEANING

All products have an approved cleaning label which you should follow to ensure durability and cleanliness of the product.

ALL WATERPROOF COVERS wipe over with hot water and soap, a mild hypochlorite solution can also be used (Milton) if desired. Vapour permeable material, towelling and spacer covers can be machine washed at 60 degrees and tumbled dry on warm only. Leave to air dry before refitting. Do not use furniture polish on the cover.

CLEANING INSTITUTIONAL VAPOUR PERMEABLE MATERIAL

Can be washed at 60 degrees Celsius or dry cleaned. Wipe over with full strength Sodium Hypochlorite or Sodium Dichloroisocyanurate and maintain surface contact for 2-5 minutes. Wipe off and rinse with clean water and dry.

DO NOT USE PHENOL BASED PRODUCTS ON ANY WATERPROOF COVER

If MRSA is suspected or confirmed please contact your Infection Control Department for advice.