

We care more



ZITZI GUPPY



Small, light
and flexible

MMS
MEDICAL



ZITZI GUPPY

Easy to balance

Balance has been the focus throughout development, both for simple testing as well as for opportunities to combine many different types of seating systems.

Multi function

Guppy is designed to fit most seating systems on the market.

Wheels

Different wheels to many needs.



QUICKIES GUPPY

SIZES.....	2
WHEELS.....	16", 20", 22" och 24"
LOAD WEIGHT.....	70 kg
TILT.....	0-32°
SEAT HEIGHT.....	46 cm
COLOR.....	Black
WEIGHT.....	14 kg
CRASH TESTED.....	Approved

“Small, light and flexible.”



Accessories Zitzi Guppy



Footrest, size 1
turn-up



Footrest, size 2
turn-up



Calf support size 1



Cushion Calf support
size 1



Bracket parkin brake
size 16"/20"/22"/24"



Parking brake



Anti tipper




Step bar




Mud flap
20"/22"/24"



Armrest



Mounting plate
X-strong size 1/2/3



Mounting system
size 1/2/3




Adapter system
Universal



Transportation kit

Technical data Zitzi Guppy



with 16" wheels (6" castorwheels)

Size	Tilt	Seat height	Weight	A	B	C	D	E
Size 1	0-32°	46 cm	14 kg	51 cm	67-73 cm	41 cm	57 cm	41 cm
Size 2	0-32°	46 cm	14 kg	56 cm	67-73 cm	41 cm	57 cm	46 cm

with 20" / 22" / 24" wheels (6" castorwheels)

Size	Tilt	Seat height	Weight	A	B	C	D	E	E with mudguard
Size 1	0-32°	46/48/50 cm	15 kg	57 cm	71-82 cm	41 cm	57 cm	43,5 cm	39,5 cm
Size 2	0-32°	46/48/50 cm	15 kg	62 cm	71-82 cm	41 cm	57 cm	48,5 cm	44,5 cm

The total length of the chassis change depending on location and wheels. C-measurements is in height with the MH for the casterwheel.