


EN

Technical data

Guppy



With 16" wheels (6" casters)

Size	Tilt	Seat height	Weight	A	B	C	D	E
1	0° - 32° -7° - 25°	46 cm	14 kg	51 cm	67-73 cm	41 cm	57 cm	41 cm
2	0° - 32° -7° - 25°	46 cm	14 kg	56 cm	67-73 cm	41 cm	57 cm	46 cm

With 20"/22"/24" wheels (6" casters)

Size	Tilt	Seat height	Weight	A	B	C	D	E	E with mud flaps
1	0° - 32° -7° - 25°	46 cm	15 kg	57 cm	71-82 cm	41 cm	57 cm	43.5 cm	39.5 cm
2	0° - 32° -7° - 25°	46 cm	15 kg	62 cm	71-82 cm	41 cm	57 cm	48.5 cm	44.5 cm

The total length of the chassis change depending on location and wheels.
C measurement is taken at the height of castor house.