



“Kids are tough on their sports chair, the Club Sport Adjustable takes it all in it’s stride! The simple seat height adjustment makes it an ideal choice for junior sports clubs, like here at Sheffield Steelers.”

► Colin Price, Chairman Steelers, Wheelchair Basketball Club

CLUB SPORT

ADJUSTABLE

The Club Sport Adjustable is the ideal solution if you are looking for a multi-sport chair with a variable seat height option.

The club Sport Adjustable is an ultra versatile chair ideal for multi-sports club or individual use. The unique adjustment system requires no tools and can be adjusted in less than 30 seconds!

► Standard Features


- ✓ Standard size option chrome moly frame (over 25 size permutations)
- ✓ Ergonomic seat
- ✓ Seat height adjustment system
- ✓ Quick release offensive wings
- ✓ Adjustable camber bar
- ✓ Adjustable height and angle footrest
- ✓ Quick release single anti-tip
- ✓ Aluminium sideguards
- ✓ Adjustable height backrest
- ✓ 5cm cushion

► Wheel Sizes

Rear wheel size: 24, 25"
Castor wheel size: 72mm

► Optional Extra's

- Click strap
- 100mm adjustable sport strap



► No tools required for adjusting the height of the chair



► Length adjustable quick release anti tip



► Adjustable height backrest no tools required



► Quick release offensive wings

If you're ready to be measured up and order your **Club Sport Adj** chair call **MMS Medical** on **021 4618000** for a free consultation.

www.mmsmedical.ie



MMS Medical, 51 East Gate Drive, Little Island, Cork
T: 021 4618000 F: 021 4618099
E: info@mmsmedical.ie

www.rgklife.com