

International
WHEELCHAIR
TENNIS OPEN



“The stabilising system on the Grand Slam is better than anything else I have come across. Combined with the custom spec of my Grand Slam I am getting the results I am after. My turning is better than ever and I can not believe a chair could help my performance as much as it has!”

► **Jamie Burdekin GB men tennis player**

grand slam

The ultimate tennis chair with a design that is great for all court surfaces and maximising performance.

The integrated camber-bar combined with a choice of fixed or adjustable footrest ensure a rigid and versatile product.

► Standard Features

- ✓ Bespoke chrome moly or titanium frame
- ✓ Integrated camber bar
- ✓ Adjustable or fixed footrest
- ✓ Length adjustable anti-tip
- ✓ Adjustable knee support system
- ✓ Adjustable tension upholstery
- ✓ Calf strap
- ✓ Foot strap
- ✓ Tool kit

► Optional Extra's

- Sideguard options
- Custom strapping
- Wheel upgrades
- Cushion options
- Upholstery upgrades
- Double anti-tip
- Taper / ergonomic seat

► Wheel Sizes


Rear wheel size 24, 25, 26", 700c
 Castor wheel sizes: 50,61,72mm



► RGK knee support system



► Titanium frame option



► Welded sideguards



► Length adjustable anti-tip

If you're ready to be measured up and order your **Grand Slam chair** call **MMS Medical** on **01 4618000** for a **free consultation**.

www.mmsmedical.ie



MMS Medical, 51 East Gate Drive, Little Island, Cork
 T: 021 4618000 F: 021 4618099
 E: info@mmsmedical.ie