

RGK

 **FRONTWHEEL**



MMS
MEDICAL




FRONTWHEEL



200x50 pneumatic tyre
with anti-flutter fork



Quick release handles to
activate and release



Robust three point
connection to frame



Retro fit to any TIGA/TIGA fx
via mounting clamps

Pushing the boundaries, see where it will take you

Developing solutions to improve your way of living is what we thrive on everyday which is why we have designed the FrontWheel.

It enables you to travel across grass, cobbles, mud, gravel and other uneven ground easily by keeping your front castor wheels off the floor at all times when you move in any direction, including backwards.

KEY FEATURES

- Exclusive to the TIGA and TIGA fx
- Three point frame connection
- Compact and easy to install in seconds
- Fits to frame with support under footrest for strength
- Castor wheels remain off the ground at all times in any direction
- Footrest width is unaffected
- No need to move feet during installation
- Locks up and down for when both in use and not in use

