

Technical Specification

Benefits

- Available in two sizes
- Tilt in space
- Fully adjustable
- Compact & easy to fold
- Lightweight
- Pelvic belt fitted as standard
- Height adjustable ergonomic pram handle
- Shopping basket fitted as standard
- Full range of accessories available
- Crash tested to ISO7176-19
- Quick release rear wheels
- Solid PU tyres

Accessories

- Leg abductor
- Lateral thoracic supports
- Occipital headrest
- Backrest extension
- Foot Straps
- Grab bar
- Contour tray
- Harness
- Sun & rain cover

Prescription form can be downloaded from our website, www.ineedawheelchair.co.uk

PART	SIZE 1		SIZE 2	
	MIN	MAX	MIN	MAX
SEAT WIDTH (pelvic supports fitted)	200mm	270mm	260mm	320mm
SEAT DEPTH	230mm	280mm	290mm	350mm
SEAT TO FLOOR HEIGHT (upright position)	-	490mm	-	520mm
BACKREST HEIGHT (from seat base)	-	510mm	-	570mm
BACKREST HEIGHT (from seat base) with back extender	-	625mm	-	720mm
THORACIC WIDTH	120mm	210mm	170mm	240mm
THORACIC HEIGHT (from seat base)	150mm	340mm	160mm	430mm
FOOTPLATE HEIGHT	130mm	280mm	125mm	330mm
FOOTPLATE DEPTH	-	160mm	-	180mm
PUSH HANDLE HEIGHT (from floor)	910mm	1250mm	925mm	1270mm
TILT IN SPACE ANGLE	15°	50°	15°	50°
OVERALL WIDTH	-	580mm	-	640mm
OVERALL DEPTH	-	790mm	-	820mm
OVERALL BACKREST HEIGHT WITH HEADREST (from floor)	-	1070mm	-	1160mm
WEIGHT	14Kg*	-	16Kg*	-
FOLDED SIZE	800x590x500mm		830x640x500mm	
*Rear Wheels Removed				

Transportation

The Marley has passed the ISO7176-19 test and is approved for motor vehicle use with a maximum user weight of 32kg.

Upholstery

The Marley comes in black/grey upholstery which is easily removable and machine washable at 40°C



 marley



Manufactured by:

RMMIS



British design and quality



The eye-catching looks of the Marley combines simplicity and effectiveness in one design.

The Marley has been developed to meet the child's individual postural needs, with a practical design for daily family living. A versatile buggy available in two sizes with a full range of accessories that complement the many standard features.

Tilt in Space

The Marley has a continuous tilt in space range from 15° to 50° to meet individual's requirements

Shopping Basket

Under seat basket for additional storage

Foldable

Quickly and easily folded in one simple operation to a compact size ideal for both transportation and storage. Removeable rear wheels

Ergonomic Pram Handle

A height adjustable handle that can be angled and easily folded to a compact position

Lateral Thoracic Support (optional extra)

For additional trunk support easily adjustable in height and width

Seat base and Pelvic Supports

Adjustable in depth and width

Foot Supports

Height adjustable and flip up to allow for easy transfer

