

BREEZY[®]

RubiX² XL™



The lightweight bariatric wheelchair

The RubiX² XL extends the user range for heavy duty users by offering seat width up to 60 cm with 170 kg user weight.



New robust castor attachment offers more rigidity and a better look (0°, 3° & 6°)



Robust reliable 4 arm cross brace



Drum brake (option)



RubiX²XL

Key Product Features:

- Heavy Duty wheelchair in a lightweight frame
- Weight capacity 170 kg
- Designed to support and offer mobility to our heavier clients by 4 arm cross brace & HD upholstery
- Adjust Seat Depth by sliding cross brace 46 - 51 cm deep
- Multiple position axle plate
- Heavy Duty Tension adjustable backrest

Features:



Axle Plate



Caster Housing



Wheel Lock



Upholstery



Tension adjustable 4 strap backrest extendable from 41 to 46 cm



Robust reliable 4 arm cross brace

Technical Specifications:

Max. User Weight:	170 kg
Chair weight	from 16.8 kgs
Overall width	52 cm, 56 cm , 60 cm
Seat width	52 - 60 cm
Seat Depth	46 - 51 cm
Seat Height:	37.5 - 51 cm
Back Height	43.5 cm (41 -46 cm adj)
Armrest	Height adjustable flip back
Hangers:	70° (lower leg length adjustable from 35 - 50 cm) Angle adjustable composite footplate (+/- 15°)
Colour:	Blue
Crash Test Approved	Yes
TGA#	145834
Warranty	Limited Lifetime