

TAiMA S-GT. The lightweight model. – 2017 version –

Product details:

- String shopping bag (up to 5 kg)
- Back rest
- Tilting aid to overcome obstacles
- Walking stick holder
- Reflectors for better visibility in the dark
- Steering dampers for enhanced convenience



Rigorous
lightweight
design:

5,3 kg*

Soft polyurethane
tyres for smoother
wheel movement

* With bag, back rest and stick holder: 5,6 kg

TAiMA rollators – Product specifications.

Sizes in centimetres (cm), weights in kilograms (kg).

TAiMA	ECO	M	S-GT	M-GT	XC Luftbereifung	XC PU-Bereifung
Colour	Metallic anthracite	Metallic brown	Metallic sepia	Metallic reseda green	Metallic red	Metallic oasis green
Height	79–95	79–95	76–86	79–95	85–95	85–95
Height when folded	81	81	76	81	88	88
Seat height	59	59	51	59	61	61
Seat width	45	45	39	45	45	45
Length	68	68	62	68	75	75
Width	63	63	58	63	67	65
Width when folded	23	23	23	23	27	25
Max. user weight	150	150	150	150	150	150
Height-adjustable pushing handles	5 x 4	5 x 4	5 x 2.5	5 x 4	5 x 2.5	5 x 2.5
Wheel size	20 x 3.5	20 x 3.5	20 x 3.5	20 x 3.5	31 x 5	31 x 5
Max. bag load	5	5	5	5	5	5
Weight	7.8	6.3	5.3 / 5.6*	5.5 / 5.8*	9.1	9.8

* Excluding bag, back rest and walking stick holder / including bag, back rest and walking stick holder

MMS Medical // T: 021 461 8000 // info@mmsmedical.ie



TAiMA –
Everything a rollator should be.



Shaping mobility.



TAiMA M-GT. The lightweight model. – 2017 version –

Product details:

- String shopping bag (up to 5 kg)
- Back rest
- Tilting aid to overcome obstacles
- Walking stick holder
- Reflectors for better visibility in the dark
- Steering dampers for enhanced convenience



Rigorous
lightweight
design:

5,5 kg*

Soft polyurethane
tyres for smoother
wheel movement

* With bag, back rest and stick holder: 5,8 kg