

Jazzy air



The Jazzy Air's advanced suspension provides unique stability in raised and lower positions. Users are able to drive safely at 3.5 mph fully elevated enabling enhanced **#socialmobility**



...A new era in powered mobility.



Jazzy POWER CHAIRS



Distributed by:



@PRIDEMOBILITYUK



Crash Tested

CORK Tel: 021 461 8000 | Fax: 021 461 8099 | DUBLIN Tel: 01 450 0582 Fax: 01 456 5889
Email: info@mmsmedical.ie | Web: www.mmsmedical.ie

Introduction

Introducing the new Jazzy Air® a new era in power mobility. The advanced stability and performance of the chassis allows the Jazzy Air to be safely elevated to 10", whilst driving at 3.5 mph.

The Jazzy Air® features Active-Trac® suspension, which provides unique stability in both the raised and lowered position. With a tight turning radius of 57.15 cm (22.5"), unique Mid-Wheel 6° Drive Design and the addition of bright LED markers combines to make the Air a smoother more enjoyable journey.

#socialmobility



Jazzy Air | Seating, Colour Options & Tech Spec

Seating Options

Choice of 16" x 16", 18" x 18" and 20" x 20" contoured high-back reclining sport seat

Colour Options



VR2 drive control with LED marker switch & Buddy button for 10" seat lift



10" Seat lift, which can be operated at 3.5 Mph



Integrated footplate and footplate design



Bright LED marker lights

Technical Specifications

Drive Wheels (cm) 25.04 cm solid tyres only	Caster Wheels (cm) Front and Rear: 15.24 cm	Max Speed (mph) 4 mph (6.4 km/h)	Ground Clearance (cm) 2.54 cm At motor bracket	Turning Radius (cm) 57.15 cm (22.5")	Overall Length (cm) 110 cm with foot rigging	Overall Width (cm) 65 cm with armrests	Seating Options (inches) 16"W - 16"D 18"W - 18"D 20"W - 20"D	Batteries 2 x 12 volt 30-35 AH Deep Cycle	Available Electronics 70A PG VR2	Weight Capacity (stone) 21 Stone 6 lbs (136 k.g)
Base Weight (kg) 82 kg	Battery Weight (kg) 11.1 kg	Warranty Lifetime limited warranty on frame. 1-year limited warranty on electronics. 18-month warranty on drive control system. 2-year limited warranty on drivetrain.		Notes <ol style="list-style-type: none"> Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. Battery weight may vary based on manufacturer. 		MMS Medical 51 East Gate Drive, Little Island, Cork. MMS Medical Unit 81 Western Parkway Business Park, Lower Ballymount Road, Dublin 12				



@PRIDEMOBILITYUK

Ensure you have the most up to date information about Pride Mobility by subscribing to our social media pages.