



MMS
MEDICAL

Single tube curved stairlifts





The benefits of a Freelift stairlift from Handicare

Handicare has an established heritage of installing stairlifts. Over time our understanding of what makes a stairlift a great stairlift has grown and allowed us to develop a range that can offer a stairlift solution for most homes.

Handicare's single tube curved bespoke stairlifts are ideally suited to staircases that turn a corner, are very steep, or are particularly narrow. Available with a range of seating systems, and optional powered features, you can specify a stairlift solution that is perfectly suited to meeting your needs today and in the future.

The unit runs along a single tubular rail which is designed to hug the edge of your staircase closely whilst leaving the rest of the staircase clear for other users.

The rail system is the same regardless of which side of the stairs it is situated, or if it travels over multiple floors only the model name will change.

If the rail follows the shorter (internal) edge of the staircase your stairlift will be a Rembrandt, if it follows the longer (external) edge it will be a Van Gogh. The Vermeer is used for installations that cover more than one storey.

Safety as standard

Handicare places the highest priority on safety and reliability. Our stairlifts are fitted with multiple safety edges as standard to prevent the stairlift from moving towards any obstruction encountered during travel.

Peace of mind

Our quality assurance systems mean that you can choose a Handicare Stairlift with confidence. All our stairlifts come with a battery back up so you know that, even in the event of a power cut, your stairlift will continue to work. Your Handicare stairlift will also have audio and visual warnings to alert you if the lift is not parked in the correct position on the rail.



Rembrandt
Inside curve

Van Gogh
Outside curve

Vermeer
Multiple floors

A seat for every user

Handicare's Freelift Stairlift range offers you the choice of different stairlift seats, each of them available in a range of colours so you can choose the style and colour most suited to your needs.

Classic

The popular contoured Classic seat has an attractive leather effect finish. Made from durable, non-slip polyurethane, it is hard wearing and easy to keep clean.

Upholstery Colour Choices



SLATE

BORDEAUX
RED


SAND



Elegance

The new leather effect Elegance seat is lightly padded and offers a harmonious combination of elegant design and sophisticated comfort.

Upholstery Colour Choices



BORDEAUX
RED

CREAM

COCOA
BROWN

Seat features



Key Switch

The key can be removed to prevent unauthorised use.



Joystick

This easy to use robust control can be positioned on the left or right arm.



Space Saving Design

The arm rests and footplate can be folded away easily, maximising the clear space on the stairs.



45° Downward Facing Seat and Angled Footplate

This seat is positioned at 45° to the wall which means that where a stairwell is particularly narrow, or a user has limited knee flexion, you can fit a Van Gogh with a 45° seat and angled footplate. This option could help avoid the inconvenience of major building adaptations, or the installation of a through the floor lift. This option is available with the Elegance and Classic seating systems.




Retractable Safety Belt

Fitted as standard for your protection and security.



Remote control

Requiring only light pressure on the call and send buttons, the lift is always within easy reach. It also allows you to park the lift at either end of the track.



Manual Swivel

The push down lever turns the seat enabling safe mounting and dismounting.



Powered Swivel*

The joystick can be used to turn the seat making it even easier to use.



Footplate Link*

Simultaneously folds up with seat.

*Optional

Handicare helping make every day life easier

At Handicare we understand that restricted mobility can cause different challenges for different people. Our aim is to provide products to help overcome as many of these challenges as we can by providing safe and reliable quality solutions worldwide.

A tailor made solution

Every Handicare curved stairlift is designed to suit the individual staircase. We will use advanced computer technology to ensure we design a rail that is the best possible fit for your staircase.

Drop nose


Handicare's tight-fitting drop nose rail can save having to fit a hinged track at the bottom of the stairs.



Automatic folding hinge

To avoid causing an obstruction at the bottom of the stairs where there is a doorway for example, an automatic hinged rail with safety system is available.

This will automatically raise and lower the bottom section of the rail thereby removing any obstruction at ground level.



Choice of rail colours

With a choice of three standard rail colours; cream, brown and light grey, you can be sure to find an option to suit your décor.



Also available from MMS Medical

Bathing solutions

Including bath lifts with a choice
of fixed or recline positions



Mobility solutions

Including award winning
walking aids

MMS

MEDICAL



Please be aware that actual colours and specifications may vary slightly from those shown. All information is provided in good faith but does not constitute a guarantee. Handicare Accessibility reserves the right to alter designs and specifications without prior notice.

MMS Medical,
51 East Gate Drive,
Little Island,
Cork

T: 021 4618000
F: 021 4618099
E: inof@mmsmedical.ie

MMS Medical,
Unit 81, Western Parkway Business Park,
Lower Ballymount Road,
Dublin 2

T: 01 4500582
F: 01 4565889
E: sales@mmsmedical.ie

www.mmsmedical.ie
www.stair-lifts.ie