



QUANTUM®

The Rehab Power Chair Solutions Company

Exceptional turning radius and modular design. The New rear wheel drive chair from Quantum Rehab.

ASPEN



Distributed by:



@QUANTUMREHABUK



Crash Tested



A cost effective rear wheel drive chair, which is a modular state-of-the-art design that provides improved obstacle climbing performance and a smoother ride.


Featuring the new TB Flex seating system, an easily adjustable, versatile design that can grow with the client and can accept Q-Logic and R-Net controls.




Colour Options



Red Corvette



Back in Black



Electric Blue

Technical Specifications

Drive Wheels (cm) 35.56 cm	Caster Wheels (cm) 22.86 cm	Max Speed (mph) 4mph	Ground Clearance (cm) 4.49 cm	Turning Radius (cm) 64.13 cm	Overall Length (cm) 107 cm	Overall Width (cm) 58 cm	Seating Options (inches) TB Flex 15" - 20"	Batteries x 2 55 AH Deep Cycle	Available Electronics VR2/RNet QLogic 2	Weight Capacity (stone) 21 Stone 6 lbs (137kg)
--	---	------------------------------------	---	--	--	--	---	--	--	---

Seat to ground height (") 17"-18.5"	Base Weight (kg) 129.48 kg	Battery Weight Each (kg) 17.21 kg
---	--	---

Warranty
 Lifetime limited warranty on frame.
 1-year limited warranty on electronics.
 18-month warranty on drive control system.
 2-year limited warranty on drivetrain.

Notes

- Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
- Battery weight may vary based on manufacturer.

MMS Medical
 51 East Gate Drive,
 Little Island,
 Cork.

MMS Medical
 Unit 81
 Western Parkway Business Park,
 Lower Ballymount Road,
 Dublin 12



@QUANTUMREHABUK