

Wheelchairs from DIETZ – Professional products for comfortable mobility.

Distributed by:

MMS
MEDICAL



»Mobility is a part of me. And I never compromise!«



NEW

ASX [01] – The adaptive wheelchair for the individual requirements of customers with geriatric and neurological conditions

The demands of the market, statutory health insurance contracts and [in Germany] the directory of aid codes are growing all the time. Accordingly, manufacturers are constantly having to modify their products and functions. Predictably, this has given rise to a virtually unmanageable array of parts and components.

The ASX[01] marked DIETZ' first bold step beyond the realms of geriatric care. At REHACARE 2014, the company thus unveiled its adaptive wheelchair for neurological care applications. The extremely adaptable ASX[01] maximizes the user's comfort and mobility, ensuring that each and every user can make optimal use of their own capabilities. Each wheelchair is manufactured to tried and trusted DIETZ quality standards on a specially designed assembly line at the new factory building in Germany.

Effective product system

DIETZ opted for a pragmatic approach and planned for a diverse spectrum of variants right from the word go. The result: A cleanly designed and well-thought-through product. Numerous adjustable settings

and configurations allow for no fewer than 22 different combinations of seat depth and seat width, for example. Combining these features with six back rest options and various accessories (leg rests, side panels, etc.) created a product concept that is truly without equal.

Efficient care

For specialist retailers providing patients with their first wheelchair, it is important to keep things simple and identify the right model with as few parameters as possible. With the ASX[01], this can be done with just two or three basic criteria (body size/weight, pelvic width and seat depth/seat width). All other parameters (front and rear seat heights, seat angle, back height, wheel camber, etc.) can be adjusted on the wheelchair itself without the need to fit extra parts.

As a result, you normally only have to visit the patient once – which saves you precious time.







The ASX[01] product concept: Every conceivable variant – in a single product.

Product details (standard configuration)

- n Aluminium frame and cross strut
- n Adjustable back rest angle: -5°/0°/5°/10°/15°
- n Adjustable upholstered back rest (Velcro)
- n Adjustable drive wheel camber: 0°/3°
- n Optional leg rest angles: 70°/80°
- n Choice of three different types of side panels (clothes guard, PU arm rests, desk & long)
- n Extremely stable front-wheel suspension
- n Ideally suited to add-on motors
- n All kinds of settings and adjustments (such as front/back seat height and back rest angle) without the need to fit extra parts

All standard configurations and optional components are listed on our interactive ASX[01] order form. Your regional DIETZ sales manager will be happy to advise you.



Atomic umbra metallic**



Petrol metallic**

** For technical reasons, printed colours may differ from the original colours.

	AS [01]
Colour	Frame matt black, side panels pearl white / petrol metallic / atomic umbra metallic
Seat depth (cm)	36 – 57 in increments of 3 48 – 57 in increments of 3
Seat width (cm)	38 – 54 in increments of 2
Seat height front (cm)	41,25 – 52,5
Seat height back (cm)	38,75 – 52,5
Back height (cm)	35 – 50 in increments of 3
Back rest angle (degrees)	adjustment range – 5° / 0° / 5° / 10° / 15°
Height of push handles (cm)	86 – 98 depending on the back height (option: extra 10 cm adjustment range)
Overall width (cm)	Seat width + 20
Width when folded (cm)	28 with wheels / 30 with wheels and drum brake
Overall height (cm)	max. 98 given a back height of 50
Overall length with leg rests (cm)	Seat depth + 45
Overall length without leg rests (cm)	Seat depth + 33
Height of arm rests (cm)	22 – 26
Width of arm rests (cm)	6
Length of arm rests (cm)	25 Desk / 35 Long
Length of leg rests (cm)	42 – 51 in increments of 1.5
Drive wheel diameter (inches)	24 / 22
Guide wheels (cm)	12,5 / 15 / 17,5 // 5 / 6 / 7
Tipping stability (degrees)	10°
Max. user weight (kg)	125 / 150
Weight (kg)	from 12,8

Wheelchairs from DIETZ – Professional products for com- fortable mobility

DIETZ' development of the professional CANEO series combines several benefits in one and highlights the company's approach to manual wheelchairs:

- ➔ Full compatibility across all components in the series
- ➔ Adjustable for all seat heights with one wheel size
- ➔ Extremely stable design thanks to innovative frame geometry
- ➔ Optimal solution for flat-rate-per-case situations and reuse

Shown: CANEO_S with drum brake



Manual wheelchairs at a glance

Our top-of-the-range model for the most exacting demands

CANEO_S

Lightweight folding wheelchair

>>> Page 9



Our oversized model: Specially for body sizes over 180 cm

CANEO_L

Lightweight folding wheelchair

>>> Page 12



Our basic model, complete with aluminium frame

CANEO_E

Lightweight folding wheelchair

>>> Page 9



Our strong model, for loads up to 200 kg

CANEO_200

Lightweight folding wheelchair for heavy loads

>>> Page 12



Our sturdy model: Outstanding value for money

CANEO_B

Standard folding wheelchair

>>> Page 10



Tokin

Lightweight folding wheelchair

>>> Page 16

Porter

Wheelchair for transportation

>>> Page 16

Serena II

Care/multifunctional wheelchair

>>> Page 17

CANEO_S

Lightweight folding wheelchair



Seat widths from 36 cm

In addition to its comprehensive standard equipment, the top-of-the-range CANEO_S model also comes with all sorts of adjustment options, enabling this wheelchair to be quickly and individually customised. With the low weight conferred by its aluminium frame, it offers unprecedented stability and conforms to the highest standards of user comfort.

- n Corrosion-resistant aluminium frame
- n Tipping stability up to 10°
- n Back angle 3°, seat angle 3°
- n Adaptable back (Velcro®)
- n Adjustable arm rests in Desk and Long versions, height-adjustable
- n Flip-back side panels
- n Adjustable lower leg length in 15 mm steps from 42 to 51 cm
- n Swing-away/detachable leg rests
- n Height-adjustable and angle-adjustable foot plates with heel strap
- n Locking brake / drum brake (depending on model)
- n Rear wheel 24", PU incl. quick-release axle
- n Covers made of high-quality plastic
- n Height-adjustable pushing handles

Accessories >>> Pages 14–15 and 18–19



The castor and adapter plate on the rear wheel enable the wheelchair to be adjusted for four seat heights without any further accessories.



Simple conversion: change from Long to Desk version in just a few seconds



Foot plates angle-adjustable

	CANEO_S					
Colour	Ivory	Ivory	Ivory	Ivory	Ivory	Ivory
Seat width (cm)	36	39	42	45	48	51
Seat depth (cm)	41 / 44	41 / 44	41 / 44	41 / 44	41 / 44	41 / 44
Seat height (cm)	42–51	42–51	42–51	42–51	42–51	42–51
Back height (cm)	40 / 43	40 / 43	40 / 43	40 / 43	40 / 43	40 / 43
Arm height (cm)	24 / 26 / 28	24 / 26 / 28	24 / 26 / 28	24 / 26 / 28	24 / 26 / 28	24 / 26 / 28
Overall dimensions W* x H x L (cm)	56 x 92 x 103	59 x 92 x 103	62 x 92 x 103	65 x 92 x 103	68 x 92 x 103	71 x 92 x 103
Overall length without leg rests (cm)	77	77	77	77	77	77
Width when folded (cm)	31	31	31	31	31	31
Max. user weight (kg)	120	120	120	120	130***	130***
Weight** (kg)	16.3	16.6	16.9	17.3	18	18.2
Item no.	210136	210139	210142	210145	210148	210151
Item no. Drum brake	211136	211139	211142	211145	211148	211151

* plus 1 cm with drum brake ** complete with leg rests and arm rests *** complete with double cross bar

CANEO_E

Lightweight folding wheelchair



Combined side panel Long/Desk

The specification of the entry level model CANEO_E – the linchpin of the CANEO series – satisfies all day-to-day care requirements. Thanks to its lightweight aluminium frame, it also offers unprecedented stability.

- n Corrosion-resistant aluminium frame
- n Tipping stability up to 10°
- n Back angle 3°
- n Seat angle 3°
- n Adjustable arm rests in Desk and Long versions
- n Flip-back side panels
- n Adjustable lower leg length in 15 mm steps from 41 to 51 cm
- n Swing-away/detachable leg rests
- n Height-adjustable foot plates with heel strap
- n Locking brake / drum brake (depending on model)
- n Rear wheel 24", PU incl. quick-release axle
- n Front wheel 7", PU
- n Covers made of high-quality plastic

Accessories >>> Pages 14–15 and 18–19



Side panels can be set to a height of 23 or 25 cm with the enclosed spacer elements.



Foot plates angle-adjustable

	CANEO_E				
Colour	Ivory	Ivory	Ivory	Ivory	Ivory
Seat width (cm)	39	42	45	48	51
Seat depth (cm)	41 / 44	41 / 44	41 / 44	41 / 44	41 / 44
Seat height (cm)	42–51	42–51	42–51	42–51	42–51
Back height (cm)	40 / 43	40 / 43	40 / 43	40 / 43	40 / 43
Arm height (cm)	23	23	23	23	23
Overall dimensions W* x H x L (cm)	59 x 92 x 103	62 x 92 x 103	65 x 92 x 103	68 x 92 x 103	71 x 92 x 103
Overall length without leg rests (cm)	77	77	77	77	77
Width when folded (cm)	31	31	31	31	31
Max. user weight (kg)	120	120	120	130***	130***
Weight** (kg)	16.4	16.7	17.1	17.8	18
Item no.	220139	220142	220145	220148	220151
Item no. Drum brake	221139	221142	221145	221148	221151

* plus 1 cm with drum brake ** complete with leg rests and arm rests *** complete with double cross bar

CANEO_B

Standard folding wheelchair



Combined side panel Long/Desk

The standard CANEO_B model with its reliable, sturdy steel frame and two seat height settings offers outstanding value for money while at the same time satisfying all basic day-to-day care requirements.

- n Tipping stability up to 10°
- n Back angle 3°, seat angle 3°
- n Adjustable arm rests in Desk and Long versions
- n Flip-back side panel
- n Adjustable lower leg length in 15 mm steps from 41 to 51 cm
- n Swing-away/detachable leg rests
- n Height-adjustable foot plates with heel strap
- n Locking brake / drum brake (depending on model)
- n Rear wheel 24", PU incl. quick-release axle
- n Front wheel 7", PU
- n Covers made of high-quality plastic

Accessories >>> Pages 14–15 and 18–19



Option: Side panels can be set to a height of 23 or 25 cm with spacer elements. (Accessories: page 19)



CardLock security system: The professional way for public and private establishments to protect their property. More information on page 17

	CANEO_B				
Colour	Anthracite	Anthracite	Anthracite	Anthracite	Anthracite
Seat width (cm)	39	42	45	48	51
Seat depth (cm)	41 / 44	41 / 44	41 / 44	41 / 44	41 / 44
Seat height (cm)	48 / 51	48 / 51	48 / 51	48 / 51	48 / 51
Back height (cm)	40	40	40	40	40
Arm height (cm)	23	23	23	23	23
Overall dimensions (W* x H x L) (cm)	59 x 90 x 110	62 x 90 x 110	65 x 90 x 110	68 x 90 x 110	71 x 90 x 110
Overall length without leg rests (cm)	81	81	81	81	81
Width when folded (cm)	31	31	31	31	31
Max. user weight (kg)	130	130	130	140***	140***
Weight** (kg)	18.4	18.5	18.8	19.5	19.8
Item no.	110039	110042	110045	110048	110051
Item no. Drum brake	110139	110142	110145	110148	110151

* plus 1 cm with drum brake ** complete with leg rests and arm rests *** complete with double cross bar

CANEO_L

Lightweight folding wheelchair



Especially for
body heights
> 180 cm

This special model complements the intelligent CANEO product line. Its equipment corresponds to that of the CANEO_S and is designed specifically for people taller than 1.80 m.

- n Corrosion-resistant aluminium frame
- n Back angle 3°, seat angle 3°
- n Adjustable castor angle
- n Adaptable back (Velcro®)
- n Adjustable arm rests in Desk and Long versions, height-adjustable
- n Flip-back side panels
- n Adjustable lower leg length in 15 mm steps from 42 to 51 cm
- n Swing-away/detachable leg rests
- n Height-adjustable and angle-adjustable foot plates with heel strap
- n Locking brake / drum brake (depending on model)
- n Rear wheel 24", PU incl. quick-release axle
- n Front wheel 7", PU
- n Covers made of high-quality plastic

Accessories >>> Pages 14–15 and 18–19



The castor and adapter plate on the rear wheel enable the wheelchair to be adjusted for four seat heights without any further accessories.

Simple conversion: switch from Long to Desk version in just a few seconds

	CANEO_L		
Colour	Ivory	Ivory	Ivory
Seat width (cm)	45	48	51
Seat depth (cm)	47 / 51	47 / 51	47 / 51
Seat height (cm)	48–57	48–57	48–57
Back height (cm)	46	46	46
Arm height (cm)	24 / 26 / 28	24 / 26 / 28	24 / 26 / 28
Overall dimensions (W* x H x L) (cm)	65 x 100 x 111	68 x 100 x 111	71 x 100 x 111
Overall length without leg rests (cm)	85	85	85
Width when folded (cm)	31	31	31
Max. user weight (kg)	150	150	150
Weight** (kg)	18.7	19	19.3
Item no.	222145	222148	222151
Item no. Drum brake	223145	223148	223151

CANEO_200

Lightweight folding wheelchair
for heavy loads



Up to 200 kg

The heavy-duty CANEO_200 model was specially designed for the obese yet is light and foldable. The low seat height permits the use of seat cushions and the side panels are extended forwards to provide support when sitting. The advanced forward and outward position of the front wheels ensures significantly greater stability as well as adequate foot room for "walking".

- n Tipping stability up to 10°
- n Back angle 3°, seat angle 3°
- n Seat height 44 and 47 cm adjustable
- n Seat depth 40 and 44 cm adjustable
- n Adjustable back (Velcro®)
- n Swing-away reinforced tubular steel side panels coupled with the side frame
- n Advanced forward and outward position of the front wheels crossbar
- n Reinforced rear wheel linkage
- n Lower leg length adjustable in 15 mm increments from 42 to 51 cm
- n Swing-away/detachable leg supports
- n Aluminium foot rest
- n Locking brake / drum brake
- n Rear wheel 24" heavy loads, PU incl. quick-release axle
- n Front wheel 7", PU
- n Reinforced covers made of high-quality plastic

Accessories >>> Pages 14–15 and 18–19

	CANEO_200	
Colour	Ivory	Ivory
Seat width (cm)	56	60
Seat depth (cm)	40 / 44	40 / 44
Seat height (cm)	44 / 47	44 / 47
Back height (cm)	43	43
Arm height (cm)	23	23
Overall dimensions W* x H x L (cm)	77 x 89 x 113	80 x 89 x 113
Overall length without leg rests (cm)	90	90
Width when folded (cm)	34	34
Max. user weight (kg)	200	200
Weight** (kg)	27	28
Item no.	151056	151060
Item no. Drum brake	151156	151160

* plus 1 cm with drum brake ** complete with leg rests and arm rests

Safe transport of wheelchair users



Situations may arise in which it is not possible to transfer an individual to a permanently fixed vehicle seat. **The CANEO_B, _E, _S and _200 models have passed the wheelchair crash test according to ISO 7176-19.** The wheelchair restraint system meets all safety standards which permits transport in a wheelchair.

Please note: Authorisation for use only applies when using the CANEO wheelchair restraint system together with the "PROTECTOR" wheelchair attachment system, associated head rests, lap belt and diagonal shoulder belt. The "PROTECTOR" attachment system and diagonal shoulder belt need to be adapted to the particular vehicle type and can be obtained directly from AMF Bruns. The AMF Bruns "FUTURESAFE" head rest installed in the vehicle may also be used.



CardLock, the security system from DIETZ

Authorities, public and private establishments naturally want to offer people with impaired mobility the greatest possible freedom of movement in their buildings. To this end, walking aids or wheelchairs are frequently provided free of charge.

However, to ensure, that – unlike their users – this equipment cannot leave the buildings, DIETZ offers a comprehensive security system. This comprises the CardLock locking system and a CANEO_B with a seat width of 45 cm. This system protects your property simply but very effectively.

Security system comprising

- n CardLock incl. a 50 cm long security chain
- n 1 card for unlocking/locking
- n 1 CANEO_B minus drum brake, seat width 45 cm

CardLock wheelchair installation kit

(supplied)

- n Wheelchair delivered with hole on right-hand pushing handle
- n You can fit the device on the chair yourself by means of a vandal-proof bolt – Pin Torx

Chain wall fitting

(not supplied)

- n The security chain is fitted to a suitable wall either by means of a vandal-proof bolt or a Ø 6 mm nail anchor
- n Wall bracket (optional)

	CardLock/CANEO_B
Colour	Anthracite
Seat width (cm)	45
Seat depth (cm)	41 / 44
Seat height (cm)	48 / 51
Back height (cm)	40
Arm height (cm)	23
Overall width (cm)	65
Overall height with CardLock system (cm)	92
Overall length (cm)	110
Max. user weight (kg)	130
Weight (kg)	18.8
Width when folded (cm)	31
Overall length without leg rests (cm)	81
Item no.	111045
Item no. Wall bracket	197024

Extensive array of accessories for the CANEO model series



Crashtest nach ISO 7176-19
Zugelassen für den Personentransport im Fahrzeug

Product description		CANEO_S	CANEO_E	CANEO_B	CANEO_L	CANEO_200	Item no.
1 Head rest, complete with clamp mounting Please note: Not suitable for use with the wheelchair restraint system!	up to 56 cm	n	n	n	n	n	400810
2 One-handed brake	left	n	n	n	n	n	190029
	right	n	n	n	n	n	190030
One-handed brake fitted ex works	left	n	n	n	n	n	190033
	right	n	n	n	n	n	190034
3 One-handed operation for use on either side (only available without the drum brake)		n	n	n	n		197022
One-handed operation	left	n	n	n	n		197017
	right	n	n	n	n		197018
Tyres 24" PU, black anodised aluminium wheel rims and hubs, spokes made of high-alloyed steel (stainless), rear wheel can be fitted on left or right sides, connecting axle can be removed easily without tools, thereby fully retaining the chair's foldability.							
4 Folding lock, black		n	n	n	n		410017
5 Pushing handles, height-adjustable, with padding (pair)	SB 39–51 cm		n	n			410023
6 Side panel, short, height-adjustable, combination Desk + Long	left		n	n			470001
	right		n	n			470000
7 Brake lever extension for locking brake incl. nut + bolt		n	n	n	n	n	198022
8 Amputee leg rest	left	n	n	n	n	n	410082
	right	n	n	n	n	n	410081
9 BW SB leg rest, 36–40 cm, horizontal and adjustable angle	left	n	n	n			197005
	right	n	n	n			197004
BW SB leg rest, SB 42–46 cm, horizontal and adjustable angle	left	n	n	n			197007
	right	n	n	n			197006
BW SB leg rest, SB 48–54 cm, horizontal and adjustable angle	left	n	n	n			197009
	right	n	n	n			197008
10 Elevating leg rest BW SH 42–53 cm	left	n	n	n	n	n	194019
	right	n	n	n	n	n	194020
11 Wheelchair restraint system with head rest See page 15	long lap belt	n	n	n		n	195007
	short lap belt	n	n	n		n	195008
Side panel raising set for the arm rest				n			194018
Adjustable back, complete with 3 eyelets and lap belt cutout for seat width	39 cm		n	n			190005
	42 cm		n	n			190006
	45 cm		n	n			190007
	48 cm		n	n			190008
	51 cm		n	n			190009

Tokin

Lightweight folding wheelchair



- n Solid steel frame
- n Side panels can be folded up
- n Short arm rests (desk)
- n Leg support, pair, can be swivelled outwards, detachable
- n Height-adjustable and swing-away foot plates with heel strap
- n Locking brake
- n Rear wheel 24", PU, incl. quick-release axle, aluminium push rim
- n Front wheel 8", PU
- n Back pocket

Porter


Transport wheelchair

Light, small, foldable



- n Made entirely from aluminium
- n Back rest can be folded out for transport
- n Back rest with bag on reverse
- n Height-adjustable side panels
- n Swing-away leg rests
- n Height-adjustable and swing-away foot plates with heel strap
- n Locking brake
- n Brake for assistant

Tokin/Porter accessories

Product description	Item no.
 Anti-tip support	192015

Accessories >>> Pages18–19

	Tokin		Porter
	140145	140150	140146
Colour	Silver metallic		Blue metallic
Seat width (cm)	45	50	45
Seat depth (cm)	40	40	40
Seat height (cm)	51	51	50
Seat angle (degrees)	–	–	2
Back height (cm)	40	40	38
Arm height (cm)	24	24	19
Lower leg length (cm)	–	–	44–48
Arm rests L x W x H (cm)	–	–	36 x 5 x 18-24
Overall dimensions L x B*x H (cm)	107 x 63.5 x 91	107 x 68.5 x 91	107 x 64 x 96
Overall height with back rest folded out (cm)	–	–	75.5
Tipping stability (degrees)	–	–	10
Width when folded (cm)	30	30	33
Max. user weight (kg)	110	110	100
Weight** (kg)	15.66	15.84	12.36
Item no.	140145	140150	140146

* plus 1 cm with drum brake ** complete with leg rests and arm rests

Serena II

Multifunctional and comfort wheelchair

Seat width:
38–53 cm



Developed specifically to meet the high demands of day-to-day care both at home and in hospital, the standard model of the Serena II wheelchair offers maximum user friendliness. Its sturdy frame, smooth functions and individually adjustable seat position make this the first choice for people who can no longer look after themselves independently.

- n Seat and back adjustable to any angle desired
- n Detachable arm rests, adjustable for height and depth
- n Swing-away, decouplable, horizontally modifiable leg rests that can also be adjusted longitudinally
- n Standard seat height 53 cm
- n Height-adjustable pushing bar
- n Locking brake
- n Drum brake (depending on model)
- n Rear wheel 24", PU incl. quick-release axle
- n Front wheel 7", PU
- n Dartex cover
- n Head rest
- n Anti-tip wheel

Accessories >>> Pages 18–19








Serena II accessories

Product description		Item no.
 Abduction wedge		411002
 Side pads with bracket	left, right	411000
 Wheelchair table	frame 1	402271
	frame 2	402299
 Set (wheelchair table and side pads)	frame 1	411025
	frame 2	411026

	Serena II	
Colour		
Seat width (cm)	38–45 / frame 1	48–53 / frame 2
Seat depth (cm)	41–50.5	41–50.5
Seat height (cm)	53	53
Seat angle (degrees)	–8 to +30	–8 to +30
Back height (cm)	53–58	53–58
Back angle (degrees)	–9 to +30	–9 to +30
Arm height (cm)	17.5–28	17.5–28
Overall dimensions (W* x H x L) (cm)	61 x 106 x 131	69 x 106 x 131
Overall length without leg rests (cm)	93	93
Max. user weight (kg)	130	130
Weight** (kg)	35.5	38.5
Item no.	390003	390004
Item no. Drum brake	390005	390006

* plus 1 cm with drum brake ** complete with leg rests and arm rests

Manual wheelchair accessories

Product description		CANEO_S	CANEO_E	CANEO_B	CANEO_L	CANEO_200	Tokin	Porter	Serena II	Item no.
 <p>Wheelchair set (with walking stick holder, seat belt, seat cushion)</p>	SB 38–40 cm	n	n	n						404062
	SB 42–43 cm	n	n	n						404063
	SB 45–46 cm	n	n	n	n		n			404064
	SB 48–51 cm	n	n	n	n		n			404065
	SB 54 cm					n				404073
 <p>Standard seat cushion incl. cover</p>	38 x 42 x 3 cm	n	n	n						403051
	41 x 42 x 3 cm	n	n	n						403052
	44 x 42 x 3 cm	n	n	n	n		n			403053
	50 x 42 x 3 cm	n	n	n	n		n			403054
	38 x 42 x 5 cm	n	n	n						403017
	41 x 42 x 5 cm	n	n	n						403018
	44 x 42 x 5 cm	n	n	n	n		n			403019
	50 x 42 x 5 cm	n	n	n	n		n			403020
	44 x 44 x 5 cm	n	n	n	n					403100
	50 x 44 x 5 cm	n	n	n	n					403101
	44 x 47 x 5 cm	n	n	n	n					403102
	50 x 47 x 5 cm	n	n	n	n		n			403103
	50 x 50 x 5 cm	n	n	n	n		n			403059
	53 x 47 x 5 cm					n				403057
	55 x 42 x 5 cm					n				195010
59 x 42 x 5 cm					n				195011	
<p>Seat cushion cover (for incontinence)</p>	38 x 42 x 5 cm	n	n	n						403033
	41 x 42 x 5 cm	n	n	n						403034
	44 x 42 x 5 cm	n	n	n	n					403035
	50 x 42 x 5 cm	n	n	n	n					403036
	50 x 44 x 5 cm	n	n	n						403060
	50 x 44 x 5 cm				n					403061
	53 x 47 x 5 cm					n				403058
<p>Seat cushion, comfort</p>	56 x 44 x 7 cm					n				845000
	60 x 44 x 7 cm					n				845001
 <p>Seat belt, press, with metal buckle</p>		n	n	n	n	n		n	402101	
 <p>Walking stick holder with high-profile cup</p>	screw-on bracket for 25 mm tubing	n	n	n	n	n	n		402243	
 <p>Wheelchair infusion holder</p>	approx. 100 cm with 2 hooks	n	n	n	n	n		n	190016	
 <p>Anti-tip wheel</p>		n	n	n	n				402584	
 <p>Shopping bag for folding wheelchair</p>	Black	n	n	n	n		n	n	400351	

Manual wheelchairs

Electric wheelchairs

Rollators

Walkers

Scooters

Electric motors

Bath & WC

Mobility & everyday aids

Respiratory aids

XXL Products

Anti-decubitus

CustomLine (XXL configurator)

Distributed by:



MMS Medical
51 East Gate Drive,
Little Island,
Cork
T: 021 461 8000
F: 021 461 8099

MMS Medical
Unit 81
Western Parkway Business Park,
Lower Ballymount Road,
Dublin 12
T: 01 450 0582
F: 01 456 5889

E: info@mmsmedical.ie
www.mmsmedical.ie

