

HSE approved supplier



Riser Recliner Chairs

HSE Contract 9232 : Provision of
Chairs, Pressure Redistribution
and Pressure Relief Aids



Feidhmeannacht na Seirbhíse Sláinte
Health Service Executive

HSE Contract 9232

Provision of Chairs, Pressure Redistribution and Pressure Relief Aids.

PRODUCTION SPECIFICATION:

Lot	Sublot	Item Description	Supplier Name	Brand Name	Supplier Code
Riser/Recliner Chair- Dual Motor	3.1	Riser Recliner Chairs – Dual Motor 40cms	MMS Medical Ltd	Pride C1 Dual Motor Riser Recliner	PRIDEC1DM
Riser/Recliner Chair- Dual Motor	3.2	Riser Recliner Chairs – Dual Motor 46cm	MMS Medical Ltd	Pride 660 Dual Motor Riser Recliner	PRIDE660DM
Riser/Recliner Chair- Dual Motor	3.3	Riser Recliner Chairs – Dual Motor 51cm	MMS Medical Ltd	Pride D90 Dual Motor Riser Recliner	PRIDED30DM

MMS Medical Contact Details

MUNSTER	<p>MMS Medical Head Office T: 021 4618000 E: sales@mmsmedical.ie</p>
LEINSTER	<p>David Delaney Product Specialist M: 086 607 8773 E: daviddelaney@mmsmedical.ie</p>
CONNAUGHT	<p>Colm Walsh Product Specialist M : 086 816 7122 E : colmwalsh@mmsmedical.ie</p>
ULSTER	<p>– DONEGAL – Colm Walsh Product Specialist M : 086 816 7122 E : colmwalsh@mmsmedical.ie</p> <p>– CAVAN / MONAGHAN – David Delaney Product Specialist M: 086 607 8773 E: daviddelaney@mmsmedical.ie</p>

Pride C1DM Dual Motor Riser Recliner

PRIDEC1DM



Technical Specifications

Seat Width	46cm / 18"
Seat Depth	44.45cm / 17.5"
Seat to Floor	44.45cm / 17.5"
Back Width	58.42cm / 23"
Seat to Top of Back	58.42cm / 23"
Weight Capacity	25 Stone / 158.7kg

FEATURES

Dual Motor Lift Riser Recline Chair

Infinite number of recline and lounging positions

Waterfall Backrest.

Back Rest & Base Separate for ease of transportation and use.

Padded Arm Rests.

Weight capacity of 25 stone (158.7kg)

Mobile Hoist Accessible

Fabric : Water Resistant Vinyl

Colour : Dark Brown.

12 month Warranty

CE Marked

No assessment Charge

No Delivery and Set Up Charge

Extending Footboard*



User Friendly Handset with LED lighting*



USB Charging Port in Handset*



User weight capacity 158.7kg/25st*



MMS
MEDICAL

For more information

☎ 021 461 8000

✉ info@mmsmedical.ie

🌐 www.mmsmedical.ie

exclusive features on the Pride Riser Recliners from MMS Medical only.

Pride 660DM Dual Motor Riser Recliner



PRIDE660DM

Technical Specifications

Seat Width
48.26 cm / 19"
Seat Depth
50.8cm / 20"
Seat to Floor
49.53cm / 19.5"
Back Width
78.74cm / 31"
Seat to Top of Back
66.cm / 26"
Weight Capacity
25 Stone / 158.7kg

FEATURES

Dual Motor Lift Riser Recline Chair

Infinite number of recline and lounging positions.

Waterfall Backrest

Back Rest & Base Separate for ease of transportation and use.

Padded Arm Rests

Weight capacity of 25 stone (158.7kg)

Mobile Hoist Accessible

Fabric : Water Resistant Vinyl

Colour : Dark Brown

12 month Warranty

CE Marked

No Assessment Charge

No Delivery & Set Up Charge

Extending
Footboard*



User Friendly
Handset with
LED lighting*



USB Charging
Port in
Handset*



User weight
capacity
158.7kg/25st*



MMS
MEDICAL

For more information

☎ 021 461 8000

✉ info@mmsmedical.ie

🌐 www.mmsmedical.ie

exclusive features on the Pride Riser Recliners from MMS Medical only.

Pride D30DM Dual Motor Riser Recliner



PRIDED30DM



Technical Specifications

Seat Width	55cm / 22"
Seat Depth	54.61cm / 21.5"
Seat to Floor	50.80cm / 20"
Back Width	78.74cm / 31"
Seat to Top of Back	66.cm / 26"
Weight Capacity	25 Stone / 158.7kg

FEATURES

Dual Motor Lift Riser Recline Chair

Infinite number of recline and lounging positions.

Waterfall Backrest.

Back Rest & Base Separate for ease of transportation and use.

Padded Arm Rests.

Weight capacity of 25 stone (158.7kg)

Mobile Hoist Accessible

Fabric : Water Resistant Vinyl

Colour : Dark Brown

12 month warranty

CE Marked.

No Assessment Charge.

No Delivery & Set Up Charge

Extending
Footboard*



User Friendly
Handset with
LED lighting*



USB Charging
Port in
Handset*



User weight
capacity
158.7kg/25st*



MMS
MEDICAL

For more information

☎ 021 461 8000

✉ info@mmsmedical.ie

🌐 www.mmsmedical.ie

exclusive features on the Pride Riser Recliners from MMS Medical only.

Chatsworth Advance Riser Recliner



The Chatsworth Advance® riser recliner from MMS Medical sets the standard as a bespoke seating solution. Available in small, medium and large sizes, this stylish chair accommodates the needs of a wide range of users looking for both comfort and postural support. The Chatsworth Advance® can be tailored to each individual's specification in terms of seat width, height, depth and armrest height to ensure comfort and support when seated and a choice of traditional scrolled or contemporary armrests to match the chair to your own style.

The Chatsworth Advance is available with a choice of single, dual or dual tilt-in-space action. The functions are easily operated via the ergonomic handset ensuring a smooth, effortless rise and recline movement. The Chatsworth Advance® can adjust its lift and tilt for the best rising position to suit each user. It is supplied as standard in maximum tilt, lowest lift position but may easily be adjusted to achieve higher flatter lift or very high virtually flat lift upon request. This can be helpful with side transfers into bed or another chair. Please speak to our seating product advisor to discuss your requirements.



Handcrafted for carefree comfort

FEATURES

- Single, dual or dual tilt-in-space action
- Standard sizes (small, medium & large) or tailored dimensions
- Choice of scrolled (supplied as standard) or contemporary armrest style
- Choice of four backrest styles - Waterfall Back, Soft/Cushioned Back, Lateral Support Back and T-Back
- Fixed seat and legrest wrap
- Choice of Castors
- MMS Medical offer a large range of upholstery options and colours for our chairs, to suit every requirement and taste. Please contact our Customer Services Department on 021 4618000 to arrange for swatches to be sent to you to ensure you are happy with your choice of upholstery. Additionally, some vinyls and fabrics may be unsuitable for your product or environment. Our helpful Customer Services Team will be happy to assist you with your choice of upholstery.



Technical Specifications

Seat Width

Small: 470mm
Medium: 520mm
Large: 520mm

Seat Depth (Single Motor)

Small: 445mm
Medium: 470mm
Large: 495mm

Seat Depth (Dual Tilt-in-Space Motor)

Small: 445mm
Medium: 470mm
Large: 495mm

Seat Height

Small: 440mm
Medium: 470mm
Large: 495mm

Arm Height (Single Motor)

Standard: 200mm
Tailored: 150-255mm

Back Height

Small: 650mm
Medium: 700mm
Large: 700mm

Safe Working Load

158kg / 25 stone



Wendover Advance Riser Recliner



Exclusively
available at

MMS
MEDICAL

The Wendover Advance® riser recliner from MMS Medical combines contemporary styling with an easy to use motor lift system. Available in small, medium and large sizes, this elegant chair accommodates the needs of a wide range of users looking for both comfort and postural support.

A selection of back options are available to suit individuals preferences: a Soft/Cushioned Back,

Waterfall Back, Lateral Support Back or T-Back. The Wendover Advance® features simple hand control functions to quickly position the chair for optimum comfort and support. A single motor for those who prefer to maintain a seat to back angle while reclining or dual tilt-in-space for the added function of an independent legrest and backrest movement.



Handcrafted for carefree comfort

FEATURES


- Choice of single or dual tilt-in-space action
- Choice of three sizes: small, medium or large
- Contemporary armrest
- Choice of four backrest styles – Soft/Cushioned Back, Waterfall Back, Lateral Support Back and T-Back
- Fixed seat and legrest wrap
- Choice of castors.
- MMS Medical offer a large range of upholstery options and colours for our chairs, to suit every requirement and taste.

Please contact our Customer Services Department on 021 4618000 to arrange for swatches to be sent to you to ensure you are happy with your choice of upholstery. Additionally, some vinyls and fabrics may be unsuitable for your product or environment. Our helpful Customer Services Team will be happy to assist you with your choice of upholstery.



Technical Specifications

Seat Width	Small: 470mm Medium: 520mm Large: 520mm
Seat Depth (Single Motor)	Small: 445mm Medium: 470mm Large: 495mm
Seat Depth (Dual Tilt-in-Space Motor)	Small: 445mm Medium: 470mm Large: 495mm
Seat Height	Small: 440mm Medium: 470mm Large: 495mm
Arm Height (Single Motor)	Small: 190mm Medium: 190mm Large: 190mm
Back Height	Small: 650mm Medium: 700mm Large: 700mm
Safe Working Load	158kg / 25 stone



Framework Purchasers

The Framework agreement will be available to (i) the HSE and (ii) entities funded under Section 38 of the Health Act 2004 (as amended) (a Section 38 Body)

For the avoidance of any doubt, the following entities are funded under Section 38 of the Health Act 2004:

- *Adelaide & Meath Hospital Incorporating the National Childrens Hospital*
- *Beaumont Hospital*
- *Brothers of Charity (Clare)*
- *Brothers of Charity (Galway)*
- *Brothers of Charity (Limerick)*
- *Brothers of Charity (Roscommon)*
- *Brother of Charity Services South East.*
- *Brothers of Charity Southern Services*
- *Cappagh National Orthopaedic Hospital*
- *Carriglea Cairde Services*
- *Central Remedial Clinic (CRC)*
- *Cheeverstown House*
- *Coombe Women's Hospital.*
- *COPE Foundation*
- *Cork Dental School & Hospital part of University College Cork*
- *Daughters of Charity Disability Support Services*
- *Dublin Dental School & Hospital*
- *Incorporated Orthopaedic Hospital of Ireland*
- *KARE*
- *Leopardstown Park Hospital*
- *Mater Child and Adolescent Mental Health Services*
- *Mater Misericordiae Hospital*
- *Mercy Foundation Hospital*
- *Muiriosa Foundation*
- *National Maternity Hospital*
- *National Rehabilitation Hospital*
- *Our Lady's Children's Hospital*
- *Our Lady's Hospice and Care Services*
- *Peamount*
- *Rotunda Hospital*
- *Royal Hospital Donnybrook*
- *Royal Victoria Eye & Ear Hospital*
- *Saint John of God Community Services Ltd*
- *Sisters of Charity – Dublin*
- *Sisters of Charity – Kilkenny*
- *South Infirmary Victoria University Hospital*
- *St. Vincent's University Hospital*
- *St. James's Hospital*
- *St. John's Hospital*
- *St. Michael's Hospital*
- *St. Michael's House*
- *Saint Patrick's Centre (Kilkenny)*
- *St. Vincent's Hospital*
- *Stewart's Care Ltd*
- *Sunbeam House Services*
- *Temple Street Children's Hospital*
- *The Children's Sunshine Home.*

OT / Client Information

Client Name:

OT:

Contact Number:

Email:

Deliver to:

Delivery Instructions

For further information or placing orders please contact

MMS Medical Contact Details

MUNSTER	MMS Medical Head Office T: 021 4618000 E: sales@mmsmedical.ie
LEINSTER	David Delaney Product Specialist M: 086 607 8773 E: daviddelaney@mmsmedical.ie
CONNAUGHT	Colm Walsh Product Specialist M : 086 816 7122 E : colmwalsh@mmsmedical.ie
ULSTER	– DONEGAL – Colm Walsh Product Specialist M : 086 816 7122 E : colmwalsh@mmsmedical.ie – CAVAN / MONAGHAN – David Delaney Product Specialist M: 086 607 8773 E: daviddelaney@mmsmedical.ie

51 Eastgate Drive
Little Island
Cork, Ireland
tel 021 461 8000
fax 021 461 8099

Unit 81
Western Parkway Business Park
Lower Ballymount Road
Dublin 12, Ireland
tel 01 450 0582
fax 01 456 5889

email info@mmsmedical.ie

www.mmsmedical.ie

MMS
MEDICAL