



QUICKIE
Live without limits.

MMS
MEDICAL

NEON²

**STYLE MEETS
PERFORMANCE.**



NEON²

STYLE MEETS PERFORMANCE.

For years Quickie manual wheelchairs have led the market with high-end Innovation & developments in technology for the active user. For the first time we're proud to roll this technology down to the mid-range...Introducing the Quickie Neon².

Bursting with personality, the lightweight and affordable Neon² provides an energy-efficient ride that you'll never want to be without. The sleek, effortlessly foldable open-frame fits easily in the smallest car – and this is a wheelchair you'll never want to leave behind!


Are you looking for individuality? Paint your personality all over the Neon² with a beautiful selection of colours for the frame, upholstery trims, wheels and forks. When you're done customising, it's not just the wheels that turn, heads will too!

Take a closer look at the new Quickie Neon² and we think you'll agree – you don't need to break the bank to enjoy a first class ride!

// FINE-TUNE YOUR BACK UPHOLSTERY
It's not just comfortable – Stand out from the crowd with a versatile palette of 5 upholstery trims.




// LIFE IS COLOURFUL
Tailor your Neon² with a selection of 32 eye-catching colours.



// HEAD-TURNING LOOK
There's also 5 distinctive styles to choose from for the castor fork and wheel.



// THE SPECIAL ONE
Accentuate your colour scheme further by matching the rim and hub of your 24" rear wheel with 3 anodized colour selections.



PERFORMANCE.

At the heart of the Neon² is a reinforced open-frame. The strength of the cross-brace and seat tubes provide our most rigid driving performance in a folding wheelchair. The unique design reduces the flex of the wheelchair when in motion - all without sacrificing the benefits of a folding frame.

Speaking of folding, the unique open-rear of the frame allows for very compact space saving dimensions. The swing-away legrests not only ensure easy transfers in and out of your wheelchair, they're removable too. Combine it with the optional half-folding back and what you have is an extremely compact package that will fit in the smallest of cars.

The versatility of the Neon² doesn't end there. The durable aluminium axle-stem allows for simple adjustments to the rear seat height and centre-of-gravity. Need to adjust the backrest angle? It's just as easy via the aluminium backrest bracket - all without the need of being without your wheelchair.





- 1 Backrest Angle**
The right back angle adjustment is important for your ergonomic posture.
- 2 Centre of Gravity**
Enables a driving characteristic from very active to passive.
- 3 Rear Seat Height**
Exactly as it should be – adjust the height to your personal need.



// PRACTICAL PORTABILITY WITH SWING AWAY FOOTRESTS
They support easy standing transfers in and out of your wheelchair and offer very compact folded dimensions for easy portability by car.





NEON²



WHEN COMPACT
DIMENSIONS AND EASY
TRANSFERS COUNT.

OPTIONS AND FEATURES OF THE NEON²

Configure your wheelchair to your individual needs. For further information on product options please check the Neon² order form.



// ADD SOME COLOUR TO YOUR UPHOLSTERY

This very comfortable and smart backrest upholstery can be tailored by 5 different colour trims.



// SEVERAL PUSH HANDLES

Choose the push handle which fits best for you. Decide between the integrated one (as shown) or the fold-down or height adjustable one.



// ANGLE ADJUSTABLE BACK

Designer backrest solution integrated into the follow of the frame. Backrest angle from 75° to 103° in 4° increments.



// ONE ARM FORK

This one arm fork, especially in combination with the coloured aluminium forks, gives your chair the special finish.



// CARBON PLATFORM FOOTPLATE

Lightweight and stylish flip-up footplate (sideways) with angle and depth adjustment.



// CARBON SIDEGUARD

Very individual and lightweight sideguard with a unique look.



// FOLD-DOWN BACKREST

This half-folding backrest folds to the front and achieves a small wheelchair size for portability.



// SURGE® LT HANDRIMS

The Gripton rubber strip on top of the ergonomic oval shaped rim offers maximum traction.

TECHNICAL DATA.

NEON² // FOLDING WHEELCHAIR

Max. user weight:	140kg	Backrest angle:	75° to 103° in 4° increments
Frame Version:	Reinforced open frame with swing-away front	Centre of gravity (COG):	20 – 120mm
Seat width:	300 – 500mm	Rear wheels:	22" / 24" / 25" / 26" (560, 610, 635, 660 mm)
Seat depth:	340 – 500mm	Camber:	0° / 2° / 4°
Seat height front:	350 – 570mm	Chair weight:	from 11kg
Seat height rear:	330 – 530mm	Frame colours:	32
Backrest height:	250 – 475mm	Crash tested:	yes

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.co.uk or alternatively is available on request in large text.



www.facebook.com/quickiewheelchairs



MMS Medical, 51 East Gate Drive, Little Island, Cork
Te: 021 4618000 Fax : 021 4618099 E: info@mmsmedical.ie

MMS Medical Unit 81 Western Parkway Business Park, Lower
Ballymount Road, Dublin 2
Tel: 01 4500582 Fax: 01 4565889 E: sales@mmsmedical.ie

www.mmsmedical.ie