

RGK

Hi·Lite


OPEN FRAME



Adjustable centre of gravity
(integrated available)



Welded castor housing with
aluminium cap and welded
castor fork



Optional folding backrest



Height adjustable footrest
and velcro calf strap

Hi·Lite



SPECIFICATIONS

Seat Width:
250-500 mm

Seat Depth:
250-500 mm

Front seat height:
380-550 mm

Rear seat height:
360-550 mm

Back height:
150-600 mm

Backrest Angle:
-10 to +20

Camber:
0-7°

Max. User Weight:
125 kg

Product Weight:
7.5-11.5 kg

Frame Colours:
25

Manoeuvre easily throughout your daily life

The Hi Lite provides premium performance with a softer ride by having a lightweight titanium open or closed frame design. The thinner tubing offers a more elegant look (25mm) and is engineered to your individual specifications.

KEY FEATURES

- Titanium frame
- Open or closed frame design
- Adjustable centre of gravity
- Choice of castor and rear wheel sizes
- Individual specifications RGKi

CLOSED FRAME

