



NURSING

# Hygiene Pulp™



General Edition

First Comes Trust [curas.com](http://curas.com)

## New for Hospitals

---

The Patent Pending C3 emptying Bag is a unique disposable urine drainage bag system designed to safely and hygienically drain a catheterized patient's drainage bag reducing risk of spills for both patient and healthcare professional.

### Curas C3 Emptying Bags (Patent Pending)

- 20460 2L -C3 Emptying Bag with SAP by Curas™ with flexi funnel
- 20465 1.5L - C3 Emptying Bag with SAP by Curas™ with small flexi funnel
- 20467 2L- C3 Emptying Bag ( Night Bag ) with SAP by Curas™ with male connector.

### Vomit Bags (Emesis) with SAP by Curas™

The Curas™ vomit ( emesis) bags offers a convenient and simple way to manage vomitus and liquid wasted in a discreet & hygienic manner. Pre-loaded with a 9g SAP by Curas™ sachet to solidify liquids and semi-liquids, preventing splashes and contains odour for easy disposable by healthcare professionals

- 03051 1.5L Vomit Bag ( Emesis) with SAP by Curas™, ml-graduated
- 03052 1.5L Vomit Bag (Emesis) ml graduated. Recommended to be used with SAP by Curas™ (ref. 08100)

### Curas Hygiene Pulp™ products

A comprehensive range of patient friendly Curas Hygiene Pulp™ is available from MMS Medical delivering best in class, cost effective solutions for healthcare professionals

- 08034 – 1L Disposable Jug with Lid
- 08030 – 4L Disposable Wash Basin
- 08051 – 850ml Disposable Urine Bottle Square
- 08000 – 2L Disposable Bedpan Liner with lid

## New for Nursing Homes

---

### Curas C3 Emptying Bags (Patent Pending)

- 20460 2L -C3 Emptying Bag with SAP by Curas™ with flexi funnel
- 20465 1.5L - C3 Emptying Bag with SAP by Curas™ with small flexi funnel
- 20467 2L- C3 Emptying Bag ( Night Bag ) with SAP by Curas™ with male connector.

### Vomit Bags (Emesis) with SAP by Curas™

The Curas™ vomit ( emesis) bags offers a convenient and simple way to manage vomitus and liquid wasted in a discreet & hygienic manner. Pre-loaded with a 9g SAP by Curas™ sachet to solidify liquids and semi-liquids, preventing splashes and contains odour for easy disposable by healthcare professionals

- 03051 1.5L Vomit Bag ( Emesis) with SAP by Curas™, ml-graduated
- 03052 1.5L Vomit Bag (Emesis) ml graduated. Recommended to be used with SAP by Curas™ (ref. 08100)

### Curas Hygiene Pulp™ products

The diverse range of the Curas Hygiene Pulp™ products from MMS Medical provides extra comfort and safety for both patients and staff for a wide variety of procedures in a nursing home environment. Switching to the Curas disposable range is cost effective and convenient solution for care homes.

- 07000 – Disposable Plastic Bedpan & Support, 2L
- 08060 – 1.5L Disposable Urine Bottle with SAP by Curas™ & lid – X-Large
- 08051 – 850ml Disposable Urine Bottle, Square
- 08100 – SAP by Curas™, fluid absorber, 6g sachet



# C3 Emptying Bags

Non-Sterile, 2L, CleanPak™

Curas C3 Emptying Bags are specifically designed for emptying urine from close system bags (CSB) and urine measurement systems. Unlike any reusable jug or plastic container, it allows the disposing of urine in a hygienic and safe manner. Curas C3 Emptying Bag is highly recommended for the demanding and complex patient situation. The C3 bags are intended to help reduce the risk of infection by providing a closed connection during the emptying of urine. The two various of Flexi-Funnel™ are designed to fit perfectly to the bottom outlet of any CSB and standard urine drainage bag. SAP by Curas™ is incorporated in C3 bags for urine crystallization absorption; thus, ensuring safer disposal of urine.

- Innovative & practical
- Ward Pak™ dispensing box
- Minimise risk of cross infection
- Safe & convenient disposal of urine
- Flexi-Funnel™ ensures excellent closed connection
- Suitable for cytotoxic waste handling



**Patent Pending**

REF	Description	Capacity (L)	Connector Type	Non-Return Valve (NRV)	SAP by Curas™	Cross- Hanger™	Spigot	Single Packaging	WardPak™ (pcs)	Carton (pcs)
20460	C3 Emptying Bag w/SAP by Curas™	2	Flexi-Funnel™	Yes	Yes	Yes	Yes	CleanPak™	25	200
20461	C3 Emptying Bag	2	Flexi-Funnel™	Yes	No	Yes	No	CleanPak™	25	300
20465	C3 Emptying Bag w/ SAP by Curas™	1.5	Small Flexi-Funnel™	Yes	Yes	Yes	Yes	CleanPak™	25	200
20466	C3 Night bag w/ SAP by Curas™, 90cm tube	2	Small Flexi-Funnel™	Yes	Yes	Yes	Yes	CleanPak™	25	250
20467	C3 Universal Emptying Bag (Night Bag) w/ SAP by Curas™	2	Male Connector with 120cm tube	Yes	Yes	Yes	No	CleanPak™	25	200

**Spigot**

Prevents accidental spillage during transportation.

**Marked Slit Cut**

For easy disposal of bag content.

**Non-Return Valve (NRV)**

Prevents reflux of urine and contaminants.

**ml-Graduation Marking**

Measuring down to 25ml and allows clear and easy reading of urine collected.

\*Patent Pending

**Patient Name & ID Box**

Writing space for patient name and identification.

**Cross-Hanger™**

Secure placement of bag on the patient's closed system bag.

**Flexi-Funnel™**

Fits all common cross valves' bottom outlet.



**Double-Sealed**

Provides extra strength and prevent leaks.




**SAP by Curas™**

Rapidly absorbs liquid and transforms into non-spill gel.

**WardPak™**

Dispensing inner box allows easy access of individual bag and saves storage space.

**C3 Emptying Bag Compatibility**



# Vomit Bag ( Emesis)

A leading Innovator in the field of fluid absorption, the range of medical consumables from our partner Curas® provides a safe and advanced solution for the convenient and simple way to manage vomitus and liquid waste (e.g. emesis) in a discreet and hygienic manner. Strong and deep bag capacity prevents overflow while printed graduation marks facilitate indication of approximate volume. The unique “Twist & Lock” feature on the rim secures the bag content and contains odour for easy disposal.

The application of the Sap by Curas™ raw material inside the vomit bags significantly improves the compliance in the field of nosocomial control. The Sap by Curas™ option crystallizes the liquid on contact for safer handling of contaminated waste by healthcare professionals. Easy to use, with patient comfort and hygiene a priority, the range of Curas® products from MMS Medical are fast becoming the products of choice for healthcare professionals in Ireland.



## FEATURES & BENEFITS

- Safe and economical alternative to using traditional emesis basins (e.g. kidney dish).
- Wide opening fits over patient’s mouth, eliminating spill-over and reducing clean-up costs.
- Practical, discreet, easy to use and conveniently manage sickness related to gastrointestinal tract, chemotherapy, motion, food poisoning, etc.
- SAP by Curas™ ( Ref: 03051) crystallizes liquid on contact for safer handling of contaminated waste.
- Each vomit Bag (Ref: 03051) contains one sachet of SAP by Curas™ for spill & odour prevention
- Hard plastic rim makes it easy to hold and handle, thus reducing risk of contamination.
- Strong and deep bag with 1.5L capacity prevents overflow.
- Printed graduation marks for indication of approximate volume.
- “Twist & lock” feature on the rim secures bag content and contains odour for easy disposal. Bag printed with simple, illustrated instructions for use.
- Bag is latex-free and PVC-free.

Ref	Description	Capacity	SAP by Curas™	Colour	Packing	Inner Bag	Carton
03051	Vomit Bag (Emesis) with SAP by Curas™	1.5L, ml graduated	Yes	White	Inner Bag	50pcs / Bag	600pcs / carton (12 inner bags)
03052	Vomit Bag (Emesis)	1.5L, ml graduated	No	White	Inner Bag	50pcs / Bag	600pcs / carton (12 inner bags)

- Helps to reduce HAIs
- Increases convenience for nursing staff
- Improves patient comfort
- Cost effective

# Hygiene Pulp™

A patient and staff friendly solution

Proper patient care is a high priority for all healthcare employees. A comprehensive range of patient friendly Hygiene Pulp™ is now available from Curas and it provides many benefits to healthcare workers and patients throughout the world.

## Disposable Urine Bottles

Curas urine bottles come in various designs and sizes meant for male and female users, in a design that allows the urinals to be placed on the floor safely without spilling over. In fact, the high-capacity 1.5L urine bottle also functions as an emptying bottle for urine bags. Nurses can easily access the urinals from the unique dispensing carton, without being compromised by a dirty utility room access. Also available is the range of rust-free stainless steel dispensing racks to hold the carton firmly on the wall, and individual urine bottle holders to hang onto the bed-frame. SAP by Curas™ is included with the range of Curas products to prevent accidental spills and improve waste handling.





The risk of cross contamination is real and significantly increased when using re-usable alternative, irrelevant of how carefully they are washed and decontaminated. By contrast, our hygiene pulp products provide improved standards of care and help to eliminate HAIs (Hospital Acquired Infections).

### Wash Basins

Curas offers two types of detergent-proof high capacity 2L, 3L and 4L bowls specially for bed-bathing. Users can choose either the disposable medical pulp-based basin with soap/detergent proof additives for up to 2 hours of use, or hold warm water for more than 4 hours. The alternative is the general bowls series that comes with our Hygiene Control Bags™ system, which is completely water resistant and with enhanced surface integrity. Both allow for clean water supply and safe disposal of the contaminated wipes etc.



## Kidney Trays/Dishes, Emesis Bowls

Curas Hygiene Pulp™ includes the all-important kidney tray, which is manufactured in accordance with Medical Device Directives – Class 1, and is vastly used in all areas of healthcare. The products have passed microbiological testing and safe for sterilization. Our unique delivery system ensures safe use of the product in both medical and surgical wards, whilst mitigating risk of infection. The unique RoomPak™ and WardPak™ systems (with barcodes) provides convenience to the user; as it delivers the product in a clean, dust free manner and allows for quick handling and allocation within the hospital setting.



Besides being single use, hygiene pulp products provide extra comfort and safety to both patients and staff. In fact, each hygiene pulp product costs less than the total cost of handling, washing and sterilizing the re-usable alternatives made from either steel or plastic. Switching to the Curas disposable pulp range provides you extra safety and convenience, at no extra cost.



## Bed Pan Liners

The disposable bedpan liners can be used together with the disposable pulp support product, with or without the Hygiene Control Bag™ and SAP by Curas™ for safe handling and disposal. The use of this unique combination is ideal for high risk or quarantined patients, as it helps to limit the spread of HAIs (Hospital Acquired Infections). Similarly, the pulp slipper pan series can also be used with Curas plastic slipper pan support at the general ward, and patients with reduced mobility.



## General Bowls

General Purpose bowls come in 1L, 2L and 3L sizes; suitable for holding liquid contents, vomitus or oral cleaning. Additionally, all disposable bowls come with a wide and firm rim that allows for easy holding, and fits perfectly on commode chairs and thus, reduces the risk of contamination. All bowls come with graduated marks embossed for clear measuring of the content, and the entire general bowl series is supplied with a range of lids for extra protection.



The diverse range means you can use our Hygiene Pulp™ products to carry out a wide variety of procedures. The range provides great versatility for nursing staff. Products are easily stored and can be effectively disposed according to clinical waste management procedures. Good insulating properties of the product mean water is kept warm and comfortable to use.



## Measuring devices

As depicted in its name, the measuring devices come with graduation marks embossed clearly on each product. The various measuring cups and jugs allow the required amount of liquid or sample content to be measured and transported. The wide edge of the mid-stream urine sampling device allows hygienic urine collection, as it fits most toilet bowls. The jugs as well as other measuring devices, support clean urine collection from urine bags, stoma bags etc. Patients using intermittent or ileostomy solutions will also benefit from these disposable options.

## Pure Fibre Deep Trays

The disposable deep trays are made of 100% pure fibre. It is water resistant, environmentally friendly and compostable. Each product has been exposed to UV rays for additional hygiene control. The sterile and non-sterile fibre trays are used in various applications, specifically at the ward level, isolation wards, ICU and in the operating rooms. It is used for holding & transporting drugs, syringes, and IV therapy sets. Furthermore, hospitals will have the benefit of making their own procedure sets by using these disposable deep trays, as they are clean and can be easily sterilized at local sterilizing facilities.



Anxiety amongst patients is a real factor in the healthcare environment and this can be reduced by using our Hygiene Pulp™ range of products.

## Digital Scales

Curas digital scales are specially designed for weighing the many Curas Hygiene Pulp™ products, and other medical items. The illuminated LCD display allows for easy reading and the smooth surface stabilizes the pulp product placement. It has multiple indicators such as; overload, low voltage, unstable and negative weight to alert the user. As an additional feature for clinical users; Curas has a unique design of digital scales which provides illustrative descriptions of key Hygiene Pulp™ products, to enable nursing and clinical staff to work in a safe and efficient manner.



## Re-usable and Plastics Support

Lightweight, but strong and very durable Plastic Pan Supports are used to improve durability and hygienic handling of pulp pan liners. It can be cleaned in automatic washer disinfectors or hand washed. To further mitigate risk of contamination, the disposable bed pans can be used along with the Hygiene Control Bag™ draped over the reusable plastic supports. Curas products are known for their high-quality finish, allowing for easy disinfection and safe re-use, subject to local protocol including the sterilization option. The range is designed for an accurate fit with corresponding Hygiene Pulp™ products and strengthen the disposable paper liner when used in situations where the patient or users body weight must be supported.



100% Recycled Material  
100% Biodegradable  
No Colour Added

Patients prefer hygienic pulp products for many reasons. It is always comforting for a patient to know that no other contamination could have taken place before using the dish, bowl or tray. Cold and rigid stainless steel alternatives are significantly less comfortable to handle and use. In comparison, the pulp material is always clean, smooth and easy to use. In addition, pulp products are biodegradable, and made of recycled paper with no added coloring.

# Hygiene Pulp™

## First Comes Trust™

Your safety is ensured by each device being produced in accordance with the key international standards MDD2007/47/EC, ISO9001:2008 and PAS29:1999 (fluid resistance test).

At Curas we say: "First Comes Trust" and this reflects the Curas approach and commitment to deliver uncompromised quality and trusted products.

The Curas Hygiene Pulp™ range is tested to ensure all products can be easily disposed of using all types of macerators for optimal environmental waste management.



MMS Medical  
51 Eastgate Drive  
Little Island  
Cork, Ireland  
tel 021 461 8000  
fax 021 461 8099

MMS Medical  
Unit 81  
Western Parkway Business Park  
Lower Ballymount Road  
Dublin 12, Ireland  
tel 01 450 0582  
fax 01 456 5889

email [info@mmsmedical.ie](mailto:info@mmsmedical.ie)

[www.mmsmedical.ie](http://www.mmsmedical.ie)

Copyright Curas Ltd

	REF	Description	Dimension L x W x H (cm)	Quantity Pcs/Bag	Quantity Pcs/Carton
Trays & Bowls	08026	Disposable 2L Bowl w/Hygiene Control Bag™, WardPak™	ø32 x 9	25	100
	08027	Disposable 3L Wash Bowl w/Hygiene Control Bag™, WardPak™	ø30 x 10	25	100
	08028	Disposable 4L Wash Bowl w/Hygiene Control Bag™, WardPak™	30 x 26 x 9	25	50
	08030	Disposable 4L Wash Basin, WardPak™	36 x 26 x 11	50	100
	08031	Disposable 1L General Purpose Bowl w/Lid, WardPak™	ø23 x 8	100	400
	08032	Disposable 3L General Purpose Bowl w/Lid, WardPak™	ø32 x 10	40	120
	08033	Disposable 2L General Purpose Bowl w/Lid, WardPak™	ø32 x 9	50	200
	08068	Disposable Kidney Tray, 700ml, RoomPak™, Cartonless	26 x 14 x 5	25	1,000
	08069	Disposable Kidney Tray, 700ml I BASIC, Cartonless	26x 14 x 5	250	1,000
	08070	Disposable Kidney Tray, 700ml, WardPak™	26 x 14 x 5	100	600
	08071	Disposable Kidney Tray, 700ml I BASIC	26 x 14 x 5	300	600
	08202	Disposable Deep Tray, S, 100% Pure Natural Fibre, 18x12.5x3.5cm (7.1x4.9x1.4in), WardPak™	18 x 13 x 4	100	1,200
	08203	Disposable Deep Tray, M, 100% Pure Natural Fibre, 20.5x12.5x4.5cm (8.1x4.9x1.8in), WardPak™	21 x 13 x 5	100	800
	08204	Disposable Deep Tray, L, 100% Pure Natural Fibre, 27x19x4.5cm (10.6x7.5x1.8in), WardPak™	27 x 19 x 5	50	500
Human Waste Collection	08000	Disposable Bedpan Liner w/Lid, 2L, w/Lid, WardPak™	36 x 29 x 10	50	150
	08001	Disposable Slipper (Fracture) Pan Liner, 1.7L, Medium (Midi), w/Lid, WardPak™	34 x 30 x 8	50	150
	08002	Disposable Slipper (Fracture) Pan Liner, 2L, Large (Maxi), WardPak™	37 x 28 x 9	50	100
	08009	Disposable Bedpan Liner, 2L, Large, Non-Pressed, WardPak™	36 x 29 x 10	40	80
	08034	Disposable Jug w/Lid, 1L, WardPak™	26 x 19 x 9	100	200
	08035	Disposable Sampling Cup w/Lid, 200ml, WardPak™	18 x 9 x 6	85	1360
	08036	Disposable Ear Irrigation Cup (Noots Tank), 300ml, WardPak™	ø9 x 15	40	160
	08040	Disposable Mid-Stream Urine Sampler, 1L, WardPak™	33 x 21 x 10	25	250
	08050	Disposable Urine Bottle, Round, 1L, w/Lid & Stick, Dispensing Carton	25 x 11 x 12	-	120
	08051	Disposable Urine Bottle, Square, 850ml, WardPak™	19 x 11 x 12	16	96
	08052	Disposable Urine Bottle, Female, 800ml	25 x 11 x 11	-	21
	08055	Disposable Urine Bottle Square, 850ml, Dispensing Carton	19 x 11 x 12	-	117
08056	Disposable Urine Bottle, Round, 1L w/SAP by Curas™, Lid & Stick, WardPak™	25 x 11 x 12	-	28	
08057	Disposable Urine Bottle, Round, 1L w/Lid & Stick, WardPak™	25 x 11 x 12	-	28	
08058	Disposable Urine Bottle, Round, 1L w/SAP by Curas™ & Lid, RoomPak™	25 x 11 x 12	5	40	
08059	Disposable Urine Bottle, 1.5L, w/Lid, X-Large, WardPak™	26 x 10 x 13	-	28	
08060	Disposable Urine Bottle, 1.5L, w/SAP by Curas™ & Lid, X-Large, WardPak™	26 x 10 x 13	-	28	
Support Products	08100	SAP by Curas™, Fluid Absorber, 6g Sachet	6 x 3	100	2400
	08101	SAP by Curas™, Fluid Absorber, 750g Bottle	8 x 9 x 23	-	12 bottles
	08102	SAP by Curas™, Fluid Absorber, 5kg bucket	28 x 19 x 21	2 bags	1 bucket
	08110	Reusable Bedpan & Support, 2L, Large (for REF 08000), Autoclavable, CleanPak™	37 x 31 x 9	1	20
	08111	Slipper/Fracture Pan & Support, 1.7L, Medium (Midi for REF 08001), Autoclavable, CleanPak™	47 x 29 x 8	1	20
	08112	Slipper/Fracture Pan & Support, 2L, Large (Maxi for REF 08002), CleanPak™	51 x 28 x 9	1	16
	08113	Bedpan & Support, 2.5L, X-Large, w/Lid, Autoclavable, CleanPak™	36 x 31 x 11	1	6
	08120	Disposable Lid for Bedpan Liner, WardPak™	29 x 17	100	1000
	08121	Disposable Lid for Midi Slipper Pan Liner, WardPak™	25 x 20	100	700
	08150	Disposable Lid for Round Urine Bottle, WardPak™	ø6 x 2	200	2000
08190	Graduated Stick for Round Urine Bottle, 1L, WardPak™	29 x 2	500	3000	
09094	Bedpan Support, Bariatric, Disposable Pulp, Extra Heavy Duty (175-200kg), CleanPak™ (for REF 08000/REF 08009)	38 x 31 x 10	1	26	
RACK AND DISPENSER	09411	Holder for Urine Bottle, Universal, Stainless Steel, Autoclavable, Ø3.5mm	16 x 14 x 20	12	36
	09420	Disposable Urine Bottle, Holder, Universal, White Corrugated Board, I BASIC, WardPak™	43 x 24 x 0.4	30	240
	11500	Dispenser Rack for Compact Urine Bottle Carton, Wall-mounted, Stainless Steel (for REF 08056/57/59/60 & 09175/76)	46 x 36 x 33	1	5
MEDICAL SCALES	05900	Digital Scale for Clinical Use (designed for Curas Hygiene Pulp™) Round Model.	ø18 x 2	1	20
	05901	Digital Scale for Clinical Use (designed for Curas Hygiene Pulp™), Rectangle Model.	20 x 17 x 2	1	20
UROLOGY	20460	C3 Emptying Bag w/SAP by Curas™ 2L		25	200
	20461	C3 Emptying Bag 2L		25	300
	20465	C3 Emptying Bag w/SAP by Curas™ 1.5L		25	200
	20466	C3 Night Bag w/SAP by Curas™, 2L, 90cm tube		25	250
	20467	C3 Universal Emptying Bag (Night Bag) w/SAP by Curas™ 2L,		25	200
MEDICAL CONSUMABLES	03051	Vomit Bag (Emesis) w/SAP by Curas™ 1.5L ml graduated		50	600
	03052	Vomit Bag (Emesis) by Curas™ 1.5L, ml graduated		50	600
	07000	Disposable Plastic Bedpan & Support, 2L, Blue, Large	37.0x30.0x10.5	1	25