

# Saoirse



*easícare*

*The Easicare Saoirse is a standard wheelchair which is lightweight, foldable and comes with a high specification. A range of optional extras allow the chair to become highly customised and the low pricing ensures excellent value for money*

**MMS**  
MEDICAL LIMITED

The Easicare range of products includes manual wheelchairs, wheelchair accessories, overhead hoist systems and mobility products. All products evolve from a user led approach to design, and are manufactured using the latest technology and components to ensure consistently high quality. All Easicare products are backed by a comprehensive manufacturer warranty and are available in Ireland, exclusively through MMS Medical. Tel. 021 461 8000

www.mmsmedical.ie



**Cork (Head Office)**  
MMS Medical Limited  
51 Eastgate Drive  
Little Island  
Cork  
Tel. 021 461 8000  
Fax. 021 461 8099

**Dublin**  
Unit 81  
Western Parkway  
Business Park  
Lwr Ballymount Road  
Dublin 12  
Tel: 01 450 0582  
Fax. 01 456 5889

**Galway**  
Unit 20  
Royal Rock  
Ballybane  
Galway  
Tel. 091 768868  
Fax. 091 768869

**Northern Ireland**  
Unit 2  
4 Ferguson Way  
Kilbegs Business Park  
Antrim  
BT41 4LZ  
Tel: 02894 46 44 45  
Fax: 02894 46 81 82

easicare

#### STANDARD SPECIFICATION

- Choice of self-propelling or transit chair
- HD version available (23st 10lbs max user)\*
- Range of four seat widths
- Choice of two seat depths
- Height adjustable armrests
- Step tube as standard
- Anti-tip Stabiliser
- Solid rear tyres
- Solid rubber castors
- Positioning Belt - Plastic Buckle
- Adjustable front fork & Axle plate
- Short wheelbase for manoeuvrability
- Amputee Setting

#### OPTIONAL ACCESSORIES

- Elevating legrests
- Seat cushion
- Tray
- Headrest
- Positioning Belts
- Adjustable Back Upholstery
- Stump Board

#### TECHNICAL SPECIFICATION

Seat Width:	36, 40, 44, 51cm
Seat Depth:	42, 46cm
Seat Height:	45-50cm adjustable
Backrest Height:	47cm
Total Width:	57, 61, 65, 72cm
Total Length:	100cm
Total Height:	94cm
Width Folded:	31cm
Armrest Height:	20-27cm adjustable
Front Castors:	8"
Rear Wheels:	24" Self-propelling 12" Transit
Total Weight:	from 13.5 Kg
Max User Weight:	125Kg / 275lbs (single X brace) 150Kg / 332lbs (double X brace)
Included in Price:	Anti-tip, Belt, Solid Tyres
Colour:	Black

\* Subject to special order

SAOIRSE SPECIFICATIONS