



Handi Turner

Weight Limits: 18 Stone

Specification: Max lifting angle 85 degrees, Pump dimensions 10" x 9" x 6"-

Movements: Lifts the client over onto one side

Options: Possum compatible

Lease: No **Hire:** Yes

The Handi turner consists of a cotton cover with a bellow inside and is easily portable due to its lightweight nature. Simply lay the cotton cover over the top cover on your bed and tuck in the excess material between the mattress and bed base, insert the hose and press the button. Turning has never been made easier, this product has been designed with both the client and carer in mind. The Handi Turner is conveniently lightweight and easily portable.

