



SLEEK STYLE AND LUXURIOUS FEATURES



Distributed by:



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®



CORK Tel: 021 461 8000 | Fax: 021 461 8099

DUBLIN Tel: 01 450 0582 Fax: 01 456 5889

Email: info@mmsmedical.ie

Web: www.mmsmedical.ie

Features and Specification


Introduction

The Apex Rapid with CTS Suspension (Comfort-Trac Suspension) delivers luxury and convenient features like a wraparound delta tiller, high level charging port, LED kerb light, storage basket and a compact easy to transport design. The Apex Rapid's sophisticated design and abundant features make it an exceptional travel mobility choice.

Disassembles into 5 pieces



CTS Front and Rear Suspension



Also Available with a:
12 Ah Battery

4x | Other optional shroud colours available.



Yellow



Black



Orange



Pink

Available Colours

The Apex Rapid comes in two optional colours:



Red



Blue

Technical Specifications

| Weight Capacity | Max Speed | Ground Clearance | Turning Radius | Overall Length | Overall Width | Front/Rear Tyres | Total Weight w/o batteries | Heaviest piece | Standard Seating | Batteries Required |
|----------------------------|---------------------|-------------------------|---------------------------|-------------------------|--------------------------|---------------------------|----------------------------|----------------------------|------------------------|---------------------|
| Stone | | | | | | | | | | |
| 21 Stone/6lbs. (137 kg) | 4 Mph (6.4 km/h) | 3.81 cm (1.5 inch's) | 113.6 cm (44.7 inch's) | 101 cm (39.7 inch's) | 52.7 cm (20.5 inch's) | Front: 20cm Rear: 20cm | 6.23 Stone 39.6k.g. | 17.8 k.g. Front Section | Folded Seat 10.4k.g | x2 18 Ah Battery |

Warranty

Two-year limited warranty on frame
One-and-a-half-year limited warranty on Controllers
One-year limited warranty Electronic Components, battery and Charger.

Notes

Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. Battery weight may vary based on manufacturer.



MMS Medical
51 East Gate Drive,
Little Island,
Cork.

MMS Medical
Unit 81
Western Parkway Business Park,
Lower Ballymount Road,
Dublin 12

CORK Tel: 021 461 8000 | Fax: 021 461 8099
DUBLIN Tel: 01 450 0582 Fax: 01 456 5889
Email: info@mmsmedical.ie
Web: www.mmsmedical.ie

