

PRODUCT DATASHEET

Slipper/Fracture Pan & Support 1.7L, Medium (Midi for REF 08001), Autoclavable, CleanPak™

| | |
|---------------|-----------------------------|
| Brand | : CURAS® |
| Category | : Nursing |
| Range | : Bedside Toileting Support |
| Product Group | : Slipper/Fracture Pans |
| REF | : 08111 |



PRODUCT DESCRIPTION

Re-usable slipper/fracture pan support made of extra strong plastic material, with antimicrobial properties. Autoclavable.

INTENDED USE

Product is used for holding human waste and other body fluids. Also used by immobile female patients for bed-side collection of urine. When used with a Disposable Slipper Pan Liner, Medium (REF 08001), the plastic support provides a perfect fit and better hygiene, strength and stability.

FEATURES & BENEFITS

- Antimicrobial technology reduces risk of cross-contamination and bacterial growth.
- Durable, heavy duty plastic gives extra strong support for individuals weighing up to 125kg.
- Re-usable and can be repetitively washed in hot water and sterilized by autoclave.
- Smooth surface is easy to clean and disinfect.
- Perfect fit and better hygiene and stability for disposable liner (REF 08001).
- Steady with flat bottom and easy to carry when full.
- Stackable for reduced storage space.
- CleanPak™, inner bag protects product from dust and contamination before reaching its point of use.
- Fits standard rack.

PRODUCT SPECIFICATIONS

| REF | Capacity | Dimensions (L x W x H) | Average Weight/Pc | Packing |
|-------|----------|------------------------|-------------------|-----------|
| 08111 | 1.7L | 47 x 29 x 7.5 cm | 301 g ± 15 g | CleanPak™ |

Composition

- Antimicrobial additive.
- 100% pure grade, non-carcinogenic polypropylene.
- Does not contain latex, phthalates, BPA (Bisphenol A) or PCB's (polychlorinated biphenyls).

Country of Origin

- United Kingdom

Sterilization

- Non-Sterile.
- Can be repetitively washed (max. 90°C) and sterilized by autoclave (max. 137°C).

Shelf Life

- 5 years

Selected Standards & Guidelines

- MDD 93/42/EEC (CE Marking, Class I), as amended by Directive 2007/47/EC.
- ISO 13485:2003
- ISO 9001:2008
- BS EN 980:2008
- BS EN 14971:2007
- EU Directive 528:2012 Biocide

Recommended Storage Condition

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

| Usage | Disposal Method |
|------------------------------------|--|
| Disposable liner + plastic support | <ol style="list-style-type: none"> 1. If regular waste, dispose of used liner into a macerator, or in accordance with the clinical waste protocol. 2. Use disinfectant spray or wipe to clean the plastic support. Autoclave when necessary. |
| Plastic support only | <ol style="list-style-type: none"> 1. Load used supports into washer disinfectant machine and follow proper operating procedures. Autoclave when necessary. |
| End-of-life of plastic support | Dispose of or incinerate in accordance with clinical waste protocol or recycle. |

Product Symbols



PACKAGING INFO

| REF | EAN-13 Barcode | CleanPak™ | Carton | Gross Weight |
|-------|----------------|------------|-----------------|--------------|
| 08111 | 5051400081118 | 1 pc / bag | 20 pcs / carton | 6.7 kg |

| Components / Material Properties | | | Recommendations on Waste Handling | | | | |
|----------------------------------|-------------------|--------------|-----------------------------------|----------|----------|---------|-------|
| Packaging | Material | Contains PVC | Incinerate | Landfill | Macerate | Recycle | Reuse |
| CleanPak™ | PE | -- | √ | (√) | | √ | |
| Outer Carton | Corrugated Carton | -- | (√) | (√) | | √ | |
| Transport Pallet | Wood | -- | | | | (√) | √ |

√ = Highly recommended

(√) = Recommended

ADDITIONAL INFO

- HS Code: 3926909790
- Documentation available upon request. Kindly contact: info@curas.com