

PRODUCT DATASHEET

Slipper/Fracture Pan & Support 2L, Large (Maxi for REF 08002), CleanPak™

Brand	: CURAS®
Category	: Nursing
Range	: Bedside Toileting Support
Product Group	: Slipper/Fracture Pans
REF	: 08112



PRODUCT DESCRIPTION

Re-usable Slipper/Fracture Pan & Support made of PVC plastic material. Washable and reusable. Packed individually in CleanPak™ (inner bag).

INTENDED USE

Product is used as a support to the Disposable Slipper (Fracture) Pan Liner, 2L, Large (Maxi) (REF 08002) for holding human waste and other body fluids. Also used by immobile female patients for bed-side collection of urine. This plastic support provides a perfect fit to the disposable pulp liner for better hygiene, strength and stability.

FEATURES & BENEFITS

- Lightweight and durable.
- Re-usable and can be repetitively washed up to 30°C.
- Smooth surface is easy to clean and disinfect.
- Perfect fit and better hygiene and stability for disposable pulp liner (REF 08002).
- Steady with flat bottom and easy to carry when full.
- Stackable for reduced storage space.
- CleanPak™, inner bag protects product from dust and contamination before reaching its point of use.
- Fits standard rack.

PRODUCT SPECIFICATIONS

REF	Capacity	Colour	Dimensions (L x W x H)	Average Weight/Pc	Packing
08112	2 L	Blue	51 x 28 x 9 cm	310 g ± 15	CleanPak™

Composition

- Made of polyvinyl chloride (PVC).
- Does not contain latex.

Country of Origin

- United Kingdom

Sterilization

- Non-Sterile.
- Can be repetitively washed up to 30°C.

Shelf Life

- 5 years



Selected Standards & Guidelines

- MDD 93/42/EEC (CE Marking, Class I), as amended by Directive 2007/47/EC.
- ISO 13485:2003
- ISO 9001:2008
- ISO 15223-1:2012
- ISO 14971:2007
- EN 980:2008

Recommended Storage Condition

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

Usage	Disposal Method
Disposable liner + plastic support	1. If regular waste, dispose of used liner into a macerator, or in accordance with the clinical waste protocol. 2. Use disinfectant spray or wipe to clean the plastic support.
Plastic support only	Load used supports into washer disinfectant machine and follow proper operating procedures.
End-of-life of plastic support	Dispose of or incinerate in accordance with clinical waste protocol or recycle.

Product Symbols



PACKAGING INFO

REF	EAN-13 Barcode	CleanPak™	Carton	Gross Weight
08112	5051400081125	1 pc / bag	16 pcs / carton	~ 6 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
Inner Bag	PE	--	√	(√)		√	
Outer Carton	Corrugated Carton	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√ = Highly recommended (√) = Recommended

ADDITIONAL INFO

- HS Code: 3926909790
- Documentation available upon request. Kindly contact: info@curas.com