

PRODUCT DATASHEET

Bedpan & Support 2.5L, X-Large, w/Lid, Autoclavable, CleanPak™

Brand	: CURAS®
Category	: Nursing
Range	: Bedside Toileting Support
Product Group	: Bedpans
REF	: 08113



PRODUCT DESCRIPTION

Re-usable bedpan & support made of extra strong plastic material, with lid and antimicrobial properties. Autoclavable.

INTENDED USE

The plastic bedpan is used for holding human waste and other body fluids. It can be used as a standalone bedpan or as a stable hygienic support for the disposable Bedpan Liner (REF 08000). It can be conveniently slid underneath the patient to support incontinence while lying in bed and also used by immobile female patients for bed-side collection of urine.

FEATURES & BENEFITS

- Antimicrobial technology reduces risk of cross-contamination and bacterial growth.
- Durable, heavy duty plastic accommodates up to 158kg – an extra strong support even for overweight individuals.
- Wide contoured rim and horseshoe aperture provide better comfort and ease of use.
- Large, removable lid helps to contain odour and prevent spillage.
- Perfect fit and better hygiene and stability for the disposable liner (REF 08000).
- Steady with flat bottom and easy to carry when full.
- Re-usable and can be repetitively washed in hot water and sterilized by autoclave.
- Smooth, flattened surface is easy to clean and disinfect.
- CleanPak™, shrink-wrapped protects product from dust and contamination before reaching its point of use.
- Fits standard rack.

PRODUCT SPECIFICATIONS

REF	Capacity	Dimensions (L x W x H)	Average Weight/Pc	CleanPak™
08113	2.5L	35.5 x 30.8 x 10.8 cm	605 g	Shrink-wrapped

Composition

- Antimicrobial additive.
- 100% pure grade, non-carcinogenic polypropylene.
- Does not contain latex, phthalates, BPA (Bisphenol A) or PCB's (polychlorinated biphenyls).

Country of Origin

- United Kingdom



Sterilization

- Non-Sterile.
- Can be repetitively washed (max. 90°C) and sterilized by autoclave (max. 135°C).

Shelf Life

- 8 years

Selected Standards & Guidelines

- MDD 93/42/EEC (CE Marking, Class I), as amended by Directive 2007/47/EC.
- ISO 13485:2003
- ISO 9001:2008
- BS EN 980:2008
- BS EN 14971:2007
- BS 2588:1992
- EU Directive 528:2012 Biocide

Recommended Storage Condition

- Store at room temperature (4 to 25°C) and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

Usage	Disposal Method
Disposable liner + plastic support	1. If regular waste, dispose of used liner into a macerator, or in accordance with the clinical waste protocol. 2. Use disinfectant spray or wipe to clean the plastic support. Autoclave when necessary.
Plastic support only	1. Load used supports into washer disinfector machine and follow proper operating procedures. Autoclave when necessary.
End-of-life of plastic support	Dispose of or incinerate in accordance with clinical waste protocol or recycle.

Product Symbols



PACKAGING INFO

REF	EAN-13 Barcode	CleanPak™	Carton	Gross Weight
08113	5051400081132	1 pc Bedpan Support + 1 pc Lid	6 pcs / carton	~5 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
CleanPak™	PE	--	√	(√)		√	
Outer Carton	Corrugated Carton	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√ = Highly recommended (√) = Recommended

ADDITIONAL INFO

- HS Code: 3926909790
- Documentation available upon request. Kindly contact: info@curas.com

PDS/08113/040117/02