



Qutee Mattress Variator

Sizes: 3' 3' 6" 4' 6" 5' (also available in 4')

Weight Limits: Single & Three Quarter - 25 stone / 158 kilos
Double & Kingsize - 40 stone / 254 kilos

Specification: Angle of Lift 90 degrees

Movements: From Laying Down Position to Sitting Up

Options: Possom Compatible, Larger Buttons on Handset.

Lease: No

Hire: Yes

The Qutee features a side operated actuator which can be fitted either side of the bed and will allow the patient to lie completely flat. The design also allows the use of two actuators for larger clients and an added benefit of increased stability. The Qutee has an incredibly smooth lifting action and is virtually silent.

