



Savu Cot/Bed

The Savu is a paediatric bed that is both attractive and versatile. With the sturdy beech wooden frame this makes a bed that is both aesthetically pleasing and safe, with no gaps. One side has a clear panel that allows the child to look around whilst the other side folds down for nursing and care, keeping the child safe yet easily accessible. Choose your frame colour and theme from our wide range.

Sizes	Width	Mattress length	Height Standard
Internal Size	30"	60"	37"
Overall Width	35"	66"	

Weight Limits: 18 Stone

Specification: Start Height 16", Lift 16"

Movements: Head up, Knee break, Hi/Lo, Trendelenburgh & Reverse Trendelenburgh

