



Timor Riser Chair Bed

Sizes:

Width	3'6" 4" 4'6"
Length	72" or 78"
Specials	2'6" 3'

Weight Limits: 25 stone (158 kilos)

Bariatric: 50 stone

Specification: Lift 16" & 19"

Start Height: from 14"

Movements: Head up, Tilt, Leg Drop

Options: High/ Low, Headboards, Hoop, Mattress, Four Locking Casters or Adjustable Static Feet

Lease: Yes

Hire: Yes

The Timor is an alternative to hoisting or sleeping in a chair thus assisting independence. Will take you from laying to sitting and then to standing with little or no effort. The riser recliner action enables the legs to go down easily whilst supporting the back of the legs and also the spine. It has a three part mattress that bends with the movement of the bed, and beech square ends come as standard.

