



Profiling Beds



Fully Adjustable Profiling Beds

The Cumulus and Cirrus, Hi-Lo electric profiling beds from Easicare provide the choice of a Homecare style or a Nursing care style design.

The adjustable height allows a range of transfer heights for the user and working heights for the carer. The four part adjustable base provides an almost infinite range of positions to ensure the optimum level of comfort, support and pressure relief.

Whichever model is chosen, a high quality, reliable and durable bed with easy to use handset controls is assured.

Two-Part Side-Rails

The design includes 2-part side-rails which allow for the user to be assisted and easier access for carers at the same time.

Choice of Mattresses

A wide choice of mattresses is also available to suit varying user needs and ensure the right mix of support and comfort is provided.

The Easicare Cumulus and Cirrus are a superb combination of functionality and quality at an affordable price.



CE Certified

Cumulus:



Cirrus:





Homecare bed



Nursing care bed

FEATURES

- Available in Two Design Styles
- Height Adjustable
- Four Section Adjustable Profiling Base
- 2 Part Side-Rails on both sides
- Adjustable Backrest
- Adjustable Knee-break
- Four Independent Lockable Castors
- Hoist Compatible
- Easy to Use Handset Controls





AVAILABLE OPTIONS

- Heavy Duty Version
Max User Weight 170Kg
- Bariatric Version
Max User Weight 250Kg
- Wide Range of Mattresses
- Extra Width Size
- Double Size
- Extra Low
- Trendelenburg option available
- Adjustable Cable Holders

FEATURES

- Available in Two Design Styles
- Height Adjustable
- Four Section Adjustable Profiling Base
- 2 Part Side-Rails on both sides
- Adjustable Backrest
- Adjustable Knee-break
- Four Lockable Castors
- Hoist Compatible
- Easy to Use Handset Controls

SPECIFICATIONS

- Bed Frame Size 200cm x 90cm
- User Weight Capacity 145kgs
- Electric adjustable;
 - Height adjustable from 30cm to 60cm
 - Back section adjustable
 - Foot section adjustable
- 4" Braked castors
- Aluminium drop down side rails.
- Powered coated steel frame

The Easicare profiling beds are available from



www.mmsmedical.ie

www.easicare.eu

All information is subject to change without notification.
Please consult your dealer with any queries you may have. E&OE.