

BENTLEY

manual tilt



Forward pivot point provides
20° of tilt with no knee rise

- Foot and hand propel in tilted position
- Rugged construction ideal for heavy duty use
- Shock absorbers for maximum comfort and protection
- Short wheel base provides excellent maneuverability
- Transit-ready option

pdg
MOBILITY

WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

WWW.PDGMOBILITY.COM

MMS
MEDICAL



Nine standard frame colors to choose from and custom colors also available. Contact us for details.




Ideal for foot propulsion



Ideal for hand propulsion



Adjustable seat height, seat depth, back angle and rear wheel



Shock absorbers are for maximum comfort and protection

Seat width, standard 14" – 22"
Seat width, custom 23" – 32"
Seat depth, adjustable 16" – 22"
Seat height, adjustable 13" – 20"
Back angle, adjustable 90° – 122°

Back height 25" or 31"
Dynamic tilt 20°
Overall width
 12" & 16" rear wheels seat width + 9.5"
 20" & 24" rear wheels seat width + 10.5"

Overall length 28.5" – 31" (no front rigging)
Weight capacity 250 lbs, custom to 450 lbs
Chair weight 39 lbs (without front rigging, wheels, anti-tippers, armrests)
Transit option available ANSI/RESNA WC-19
Warranty PDG's exclusive TGIF warranty (Thank Goodness It Fits)
 Lifetime warranty on the frame

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular needs of the consumer and the limitations of the product.