

STELLAR

45° of tilt manual



- Ideal for hand propulsion with accessible rear wheels
- Knees remain low in tilt—wheel under desks, into vans and access other environments with ease
- Short wheelbase provides excellent maneuverability
- Features adjustable seat height, seat depth and back angle
- Available in Heavy Duty (HD) extra-wide sizes
- Transit-ready option
- ISO tested



Stellar HD

pdg
MOBILITY

WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

WWW.PDGMOBILITY.COM

MMS
MEDICAL



Nine standard frame colors to choose from and custom colors also available. Contact us for details.




Multi-adjustable



Adjustable push bar option



Easily adjustable axle position



Actuator adjusts for optimum size



HD option

Seat width (standard) 14" – 22"
 Seat width (HD) 22" – 32"
 Seat depth (adjustable) 16" – 22"
 Seat height (adjustable) 14" – 20"
 Back angle (adjustable) 90° – 120°
 Back height 18" – 26"
 Dynamic tilt 45°

Overall width
 12" & 16" rear wheels seat width + 9.5"
 20" & 24" rear wheels seat width + 10.5"

Overall length 30" – 37" (without front rigging)
 Weight capacity 250 lbs, HD to 450 lbs
 Chair weight 40 lbs (standard); 50 lbs (HD)
 Transit option available ANSI/RESNA WC-19 (Standard frame only)
 Warranty PDG's exclusive TGIF warranty (Thank Goodness It Fits)
 Lifetime warranty on the frame.

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular needs of the consumer and the limitations of the product. US Patent # 6 447064 B1

MMS Medical
 51 East Gate Drive,
 Little island,
 Cork

Tel: 021 4618000
 Fax: 021 4618099
 Email: info@mmsmedical.ie

MMS Medical
 Unit 81 Western Parkway Business park
 Lower Ballymount Road,
 Dublin 2

Tel: 01 4500582
 Fax: 01 4565889
 Email: sales@mmsmedical.ie

www.mmsmedical.ie