

FUZE JR



FUZE JR



Features

50° of manual tilt	Constant center of gravity
Adjustable axle and seat height	Short wheelbase for maneuverability
Narrow overall width	Transit tie-down
Low seat-to-floor height	



MMS Medical
51 Eastgate Drive,
Little Island,
Cork.

Tel : 021 4618000
Email : info@mmsmedical.ie
www.mmsmedical.ie

www.pdgmobility.com

The Fuze JR is a Tilt-in-Space wheelchair designed for children. With an adjustable the wheelchair frame that can be adapted to the growth of the occupant. It provides a tilt range of 50° which is virtually maintenance-free due to its unique single-point centre-of-gravity pivot. Its effortless positioning performance facilitates pressure relief, feeding and respiratory function.


A Adjustable Seat Height

13" to 20" of front seat height to optimize occupant comfort and mobility requirements.



B Adjustable Axle and Seat Position

Adjust the center of gravity location to optimize mobility requirements.




C Adjustable Seat Depth and Back Rest

A large selection of seat depth options and up to 36 degrees of back angle adjustment for general positioning and comfort.



D Adjustable Chair Width

On-board adjustability of the chair width needs to accommodate occupants growing needs.



SEAT WIDTH
12"–18"



SEAT DEPTH
13"–18"



FRONT SEAT HEIGHT
13"–20"



BACK ANGLE
84°–120°



DYNAMIC TILT
Up to 50°



MAXIMUM CAPACITY
250 lbs/113 kg



CHAIR WEIGHT
31 lbs/14 kg



OVERALL WIDTH
+8.5"

