

Iris

The Superior Standard for
Wheelchair Tilt Technology

IRIS

The best just got better. The IRIS's patented rotation-in-space technology allows for a small footprint on a lighter, technically improved, sleek design.



Shortest possible wheelbase without sacrificing stability



Lightweight

The new Quickie IRIS is one of the lightest positioning wheelchairs on the market. With more than 2.26kg of reduced weight, the IRIS is even easier to transport and maneuver.

Streamlined design

The IRIS look and function have been enhanced with our new rocker design, black hardware and welded castors.

Two tilt range options

The IRIS provides a 55° range of rotation with the standard rocker arm or a 40° range with the shorter, lightweight rocker arm.

A new foot release tilt actuator

Eliminates cables. Also available, the existing hand triggers feature a special new Cable Management System that clips the cables to the frame to limit interference and opens space for options that mount on the backrest.




Low seat-to-floor height

With our new castor design, low castor fork and small castor wheel options, the Quickie® IRIS can achieve a seat-to-floor height as low as 32 cm. This provides easier access to tables, desks and counters as well as enhanced floor access for foot propulsion.



The New IRIS
is JAY3 compatible

Short, Small, light-easy to tilt.
Unsurpassed maneuverability.



Technical data sheet

16 individual frame colours	Back height: Angle adj. stroller handle	40 - 51cm (15.7 - 20'')
Weight of heaviest part	from 16kg	Back options: Stroller handle back post (Low, medium, tall)
Max. user weight - Standard 113kg (17.7 stone)	Hanger options: 60°, 70°, 80°, 90° (Swing-away and Heavy duty) ELR, ALR	
- Heavy duty 158kg (24.8 stone)	Transit tested to ISO 7176 Part 19 and passed	Yes

For further information on the full specification, options and accessories, please refer to the order form.
All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.
If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.



MMS Medical,
51 East Gate Drive,
Little Island,
Cork

MMS Medical
Unit 81,
Western Parkway Business Park,
Lower Ballymount Road,
Dublin 2



PB_IRIS_EN_GB_REV1

Tel: 021 4618000
Fax: 021 4618099
Email: info@mmsmedical.ie

Tel: 01 4500582
Fax: 01 4565889
Email: sales@mmsmedical.ie

www.mmsmedical.ie