



**RGK**

# TIGA



TIGA folding backrest with activation cable and upholstery guides



Low ride TIGA castor fork, black anodised



Neoprene frame protectors with calf strap location guide



Integrated camber bar option

# TIGA



## SPECIFICATIONS

**Seat Width:**  
250-500 mm

**Seat Depth:**  
250-500 mm

**Front seat height:**  
380-550 mm

**Rear seat height:**  
360-550 mm

**Back height:**  
150-600 mm

**Backrest Angle:**  
-10 to +20

**Camber:**  
0, 1, 2, 3, 5°

**Max. User Weight:**  
125 kg

**Product Weight:**  
6.5-11.5 kg

**Frame Colours:**  
25

**Life in the safe lane just got faster!**

**Ideal for all aspects of busy daily life.**

The lightweight aluminium open or closed frame of the TIGA with its Airtech upholstery and forged design castor arm, and made to your individual specifications using RGKi is available with a choice of anodised accent colours.

## KEY FEATURES

- Aluminium 7020 frame
- Open or closed frame design
- Forged design castor arm
- Airtech upholstery system
- Choice of anodised accent colours
- Individual specifications with RGKi

