



AERO X



# THE TILITE AERO X. EXTENDING THE DEFINITION OF THE ULTIMATE RIDE™.

## TILITE INNOVATION.

- Exclusive **Rigid Folding Mechanism “RFM”** increases performance by reducing flex in the side frames.
- **NEW Easy Fold RFM** makes folding the Aero X a breeze.
- **Tru-Fit System:**
  - Reduces the number of holes in the frame, improving durability.
  - Allows for easy adjustment and reduces set-up time.
- **Adjustable Camber Mount** incorporates TiLite’s Tru-Fit System with up to 5” of Rear Seat Height adjustment.
- **Bullet** caster housing adjusts to accommodate the greatest range of users.
- **Performance 5-Spoke Soft Roll Wheel** weighs just 1.6 ounces, 40% lighter than the original 3-Spoke Wheel.

## TECHNICAL SPECIFICATIONS.

Frame Type	Aluminum Folding – Full Adjustability
Frame Styles	Fixed Front End Swing Away (Swing In or Out Removable) <ul style="list-style-type: none"><li>• Hemi Height Frame</li></ul> Power Adaptable Frame
Total Width Folded*	10.75”
Total Chair Width*	23”
Weight Limit	250 lbs. and 350 lbs. HD Edition
Seat Widths	12” - 22”
Seat Depths	12” - 20”
Front Seat Heights	13.5” - 22”
Rear Seat Heights	13” - 20”
Cambers	0°, 2°, 4°, 6° or 8°
Front Frame Angles	Fixed Front - 70°, 80° or 85° Swing Away - 70°, 80°, 85° or 90°
Front Wheels	24 Options
Rear Wheels	61 Options
Handrims	12 Options
Center of Gravity Adjustment	6”
Transport Weight*	Fixed Front - 16 lbs. Swing Away - 18.4 lbs.
Convenience	Transit Tie-Down Option Folding Angle Adjustable Seat Back Depth Adjustable Seat Back
Aesthetics	19 Frame Colors 7 Upholstery Colors 8 Spinerger Wheel Colors
HCPCS Code	K0005
Warranty	Lifetime on Frame; Limited Lifetime on Components

\* Weights and widths based on 16” x 16” chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24” rear wheels weighs 7.6 lbs.

## CUSTOM CONFIGURATION THE TILITE WAY.

Custom configuration is when you are able to select options and components on a chair to meet your needs. While other companies offer some custom configuration, TiLite offers the widest range of options from which to choose. More choices mean more ability to truly create a custom configuration of options that best meets your needs. Choice matters and that is why TiLite offers more options than any other manufacturer! The choice is yours.



TiLite | 800.545.2266 | [tilite.com](http://tilite.com)  
© 2011 TiSport, LLC.

All TiLite wheel chairs are shown with available options. CF\_AEROX\_0911RevC