16 BACK SUPPORTS

# **ACTA-EMBRACE™**

# **ACTA-EMBRACE™**



Also available as a custom back See page 58 for details.



E2613 / E2614 POSITIONING BACK, POSTERIOR



**FRONT VIEW** 



SHELL VIEW

**ACTA-EMBRACE™ LTS** 



E2620 / E2621
POSITIONING BACK, PLANAR WITH LATERAL SUPPORTS



### **MOLDABLE SUPPORTS**



#### MOLDABLE STAY POCKETS

# ACTA-EMBRACE™ LTS LATERAL TRUNK SUPPORTS

See page 7 for details.



# LATERAL TRUNK SUPPORT OPTIONS

See page 7 for details.

4" D x 3" L (LAT-4D3L-COMFORT-TEK or LAT-4D3L-STRETCH-AIR)

**5" D x 4" L** (LAT-5D4L-COMFORT-TEK or LAT-5D4L-STRETCH-AIR)

**6" D x 4" L** (LAT-6D4L-COMFORT-TEK or LAT-6D4L-STRETCH-AIR)

**7" D x 5" L** (LAT-7D5L-COMFORT-TEK or LAT-7D5L-STRETCH-AIR)

Lateral trunk support option can be selected by adding the (part number) to the Acta-Embrace™ LTS part number.

# ACTA-EMBRACE™ HARDWARE OPTIONS See page 4-5 for details.

Compass® 4 Quick Release 2 Point (Q-COMPASS)

Compass® 4 Fixed 2 Point (NQ-COMPASS)

Compass® 4 Quick Release 4 Point (Q-COMPASS4)

Compass® 4 Fixed 4 Point (NQ-COMPASS4)

Compass® Power Mount Bracket See page 5

Hardware option can be selected by adding the (part number) to the base back support number  $\,$ 

# **COVER OPTIONS** See page 2 for details.

Comfort-Tek™ (COMFORT-TEK)

Stretch-Air™ (STRETCH-AIR)

Cover option can be selected by adding the (part number) to the base back support number

#### WEIGHT CAPACITY

10" - 18" Wide	250 lbs
20" Wide	300 lbs
22" Wide	350 lbs
24" Wide	400 lbs

## CLINICIAN TIP

Ideal for a complex user that needs a truly custom fit

### **ACTA-EMBRACE™ FEATURES**

- Integrated trunk wings attach to a planar aluminum back support shell are adjustable for angle, height, and up to a total of 3" of width for standard wings and 1.5" for pediatric wings to create a custom supportive shape.
- Moldable stays included for positioning, support, pressure redistribution and on-the-spot customization for a variety of spinal curvatures and changes in conditions.
- Moldable stays on the backside of the two layer foam pad can be molded to support lumbar, scoliosis and kyphosis needs.
- Moldable stays are offset from center of the back, removing pressure points from the spine.
- · Comes with integrated chest harness strap guides.
- Durable aluminum shell comes with pre-drilled 2" x 1" head support hole pattern (BodiLink™ Head Support page 20).
- Compatible with lateral trunk support hardware. A total of 11" of pad width adjustment (5.5" per side) & over 5" of height adjustment. Lateral trunk supports (E0956) are optional & available on page 22.

#### STANDARD SIZES

	10" W	12" W	14" W	16" W	18" W	20" W	22" W	24" W
22" L				16 x 22	18 x 22	20 x 22	22 x 22	24 x 22
20" L			14 x 20	16 x 20	18 x 20	20 x 20	22 x 20	24 x 20
18" L		12 x 18	14 x 18	16 x 18	18 x 18	20 x 18	22 x 18	24 x 18
16" L		12 x 16	14 x 16	16 x 16	18 x 16	20 x 16	22 x 16	
14" L	10 x 14	: 14 // 11		10 / 11	: 10 / 11	20 x 14		
	10 x 12							
10" L	10 x 10	12 x 10	14 x 10	16 x 10				

\*Length (L) refers to the actual size dimension of the support from bottom to top edge.

# ACTA-EMBRACE™ PART NUMBER

AE-BSXXWYYL

# ACTA-EMBRACE™ LTS PART NUMBER

**AELTS-BSXXWYYL** 

XX & YY to be replaced / XX = width & YY = length in inches

<sup>\*24&</sup>quot; Wide back supports require 4 Point Hardware or Power Mount Hardware.

<sup>\*</sup>Power Mount Hardware not available on 10"W or 10"L back supports.

<sup>\*4</sup> Point Hardware not available on 10" Lif using head support.