


TECHNICAL DATA DELFI PRO

ADJUSTABLE BACK

All measurements are including lining



	Size 0	Size 1	Size 2	Size 3-4
Back profile length:	33	36	45	52 cm

Back dimension (cm)

	B x E
Back Size 0	34 x 21 cm
Back Size 1	34 x 21 cm
Back Size 2	39 x 26 cm
Back Size 3	45 x 29 cm

Lateral support dimension (cm)

	C x D
Lat. support Size 0	16 x 12 cm
Lat. support Size 1	20 x 13 cm
Lat. support Size 2	23 x 15 cm
Lat. support Size 3	25 x 18 cm

LATERAL SUPPORT ADJUSTMENT (A) cm incl. lining

	Back 0	Back 1	Back 2	Back 3
Lat. support size 0	17 - 26	17 - 26	19 - 28	21 - 30*
Lat. support size 1	17 - 26	17 - 26	19 - 28	21 - 30*
Lat. support size 2	19 - 28	19 - 28	21 - 30	23 - 30*
Lat. support size 3	19 - 28*	19 - 28*	21 - 30*	23 - 32

Grey boxes are standard configuration

* Not optimal height adjustment.