

QUANTUM® **KOZMO** PORTABLE



The Kozmo combines manoeuvrability with a compact design and an impressive tight turning circle of 64.5cm, making it the perfect choice for confined spaces. Packed with great features such as adjustable seating, the Kozmo can easily be dismantled for transport or storage.

UNIQUE TECHNOLOGY

to the industry.

*Small Compact Frame | Easy to Dismantle
for Transport | Smaller Seat Sizes*

MMS
MEDICAL

@Quantumrehabuk



Version 3.0 | September 2018



 **QUANTUM**

MMS Medical, 51 Eastgate Drive, Little Island, Cork. Tel : 021 4618000 Tel: 01 4500582
| Web: www.mmsmedical.ie | Email: info@mmsmedical.ie






Disassembles into 4 pieces

Kozmo | Colour Options



Ruby Red



Pearl White



Rose Quartz



Citrine Yellow



Amber Orange



Sapphire Blue

Technical Specifications

Drive Wheels (cm) 22.86 cm (9") solid	Caster Wheels (cm) 15.24 cm (6") solid	Max Speed (mph) 3.6 Mph (5.79 kph)	Ground Clearance (cm) Anti-tip: 4.06 cm (1.6")	Turning Radius (cm) Inc Footplatform 64.5 cm (25.4")	Overall Length (cm) Inc Foot Rigging 83.82 cm (33.09")	Overall Width (cm) 57.4 cm (22.6")	Seating Options 12"W x 12"D to 18"W x 18"D	Batteries x 2 18 AH Deep Cycle	Available Electronics GC3 Control System	Weight Capacity (stone) 17 Stone 8lbs 113 k.g.
Seat to ground height (") 15.75"-16.75"	Base Weight (kg) 30.4 kg	Battery Pack Weight (kg) 14.06 kg	Warranty <ol style="list-style-type: none"> Lifetime limited warranty on frame. 1-year limited warranty on electronics. 18-month warranty on drive control system. 2-year limited warranty on drivetrain. 				Notes <ol style="list-style-type: none"> Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections. Battery weight may vary based on manufacturer. 			

Get Social with Quantum



FOLLOW THE LIVES OF OUR QUANTUM STARS, AND GET THEIR PERSPECTIVE OF THE Q6 EDGE 2.0

www.quantumrehab.co.uk/quantum-stars



@QUANTUMREHABUK

Ensure you have the most up to date information about Quantum Rehab by subscribing to our social media pages

#ileveluk

#myquantumproduct

#socialmobility

Further information can be obtained through our web site, or prescription forms. All information is subject to change without notification. Please consult Quantum with any queries you may have.



Non-Crash Tested

Tel: 021 461800 or 01 4500582 | Email: info@mmsmedical.ie | Web: www.mmsmedical.ie