

Active Tray

Modular
Activity
Tray

Modular • Adjustable • Robust
Folds Down for On-Board Storage when not in use

The Modular Tabletop that Mounts at Midline

The Active Tray is a center positioned activity tray. The stability of the tray is an enormous benefit that comes from using the center mounted socket as the base, as compared to the more common side armrest mounted trays. This new accessory can be removed and conveniently stored in the side socket by the patient when not in use. The tray has multiple adjustments at almost any angle due to the ball and socket design.



Works with the Center Drive System or
independently with Manual Wheelchairs

Easy Installation • Multiple Socket Locations

MMS
MEDICAL

MMS Medical, 51 East Gate Drive, Little Island. Cork

T : 021 4618000 E : info@mmsmedical.ie W : www.mmsmedical.ie

W ActiveControls.com

**ACTIVE
CONTROLS**