




Adjustable back

For improved postural control and reduced risk of deformity while increasing potential for activity, the adjustable back allows you to correct and maintain an optimal position based on individual needs and status.




Adjustable hip/back angle & tilt

A stable base of support is needed to assist the user in creating a functional seated position. A dynamic system is also necessary in order to allow for freedom of movement, variation, and as a means of reducing the effects of extension spasticity.



Pressure reducing lining

The material used for Zitzi Clozitts lining reduces pressure thereby improving pressure distribution over an increased area. It is also specially developed for hygienic purposes as it repels both moisture- and soil.



Adjustable leg parts

The thighs supports can be adjusted exactly to suit seat depth and abduction, thus creating an optimal seated position which can “grow” with the user. It is even possible to create an open hip angle if needed.

The picture is showing
A Clozitt size 1-2 with divided seat part
and divided footrests with calf support

Headrest

With three adjustable angles you are provided with an endless number of positions for an optimal fit

Mounting bracket

To ensure the proper placement of positioning harnesses such as Sternum Support or the H-Strap

Adjustable back

For improved postural control

Driving handles

Driving bow also available

Zitzi Clozitt

Strong stainless construction

Ball rod adapter

Universal adjustment possibilities

Step bar

Easy to transport

Zitzi Clozitt Quickies

A combined hygiene- and toilet seat. Perfect for bathing in a pool! Durable construction and versatile adjustment possibilities provide continuous support and excellent positioning.

Sizes 4 st
Editions 6 st
Max user weight 130 kg
Tilt 30°

Total width 48–53 cm
Seat part, standard 26,5–41,5 cm
Divided seat part 33–47 cm
Seat height 28–80 cm



Zitzi Clozitt is the undisputed leader in flexible seating systems. Along with seamless adjustments in the seat and back, the back angle is equipped with a hydraulic spring which is helpful to those with extension spasticity.


On top of the numerous adaptation possibilities already available with the product, custom made modifications can also be requested. For advice, don't hesitate to contact our trained professionals who will be happy to assist you with ideas.



Headrest



Headrest (long)



Driving bow



Upper mounting bracket



Adjustable back



Urine protection



Armrest



Table



Pot



Footrest



Divided footrest, ball joint



Calf support (and cushion)



Raising stay



Driving handles



Anti-tipper



Sternum support



H-strap



Hip belt



Ankle holder



Bracelet



Deepening stays



Wheels



J-brackets



© Anatomic SITT AB