



CARLO Classic Alu

The ideal tool for safe and effortless lifting, repositioning and transferring of your residents who are limited in their mobility

CARLO Classic Aluminium sling hoist

offers maximum safety and comfort ...



The ideal relief in daily care ...

The CARLO Classic Aluminium is characterised by an extremely light but nonetheless very stable aluminium construction and offers maximum safety and comfort to your residents and your nursing staff.

The lifting arm with a 4-point sling bow can be lowered to the floor by means of a hand control which ensures that a person who has fallen to the floor can be safely and easily lifted by just one staff. The carrier frame can be spread electrically and is thus very comfortable to operate. Due to the low height of the carrier frame, the CARLO Classic Alu can easily pass underneath both bathtubs and care beds.

The stable sling bow is suitable for all loop slings, is provided with a safety locking mechanism to attach the slings and is characterised by a patent-pending adjustment option to reduce rotation.

Our mobile Carlo Classic Aluminium sling hoist supports your nursing staff in their everyday lifting and transfer routines in an optimal and reliable manner and allows for ergonomically safe and efficient work at all times, without putting any strain on the musculoskeletal system when lifting and repositioning the resident manually.



Benefits for your nursing staff and your residents...

- ❑ Light construction, as the carrier frame and the lifting column are made of 100% aluminium
- ❑ Extra light and easy to manoeuvre
- ❑ High-quality lockable STEINCO double castors
- ❑ Electrical height adjustment and electrically spreadable carrier frame
- ❑ No unsettling pendular or rotational movements during lifting and transferring, with patent-pending adjustment option to reduce rotation
- ❑ Hand control with display of the battery charge status
- ❑ Working in an ergonomic and health-orientated manner, as no manual lifting is required
- ❑ 4-point sling bow with safety locking mechanism
- ❑ Battery and charger compatible with all BEKA Hospitec hoists
- ❑ Large choice of comfortable loops slings of different sizes
- ❑ Also available with patient scale



Technical data of the CARLO Classic Aluminium

Carrier frame (length/height)	1.290 mm / 125 mm
Inner width (closed/spread)	540 mm / 1.040 mm
Outer width (closed/spread)	700 mm / 1.180 mm
Total width of the lift (carrier frame closed)	700 mm
Lifting range	1.400 mm
Max. patient weight	185 kg
SWL	185 kg
Drive	24V LINAK drive
Power consumption	1 ampere

Accessories/Options

Digital patient scale
Transfer slings with loops (available in different sizes, standard sizes S-XXL)
Bathing slings with loops (available in different sizes, standard sizes S-XXL)
Toilet slings with loops (available in different sizes, standard sizes M-XL)
Spare battery 24V
External wall charging unit



MMS Medical
51 East Gate Drive,
Little Island,
Cork

T: 021 4618000
F: 021 4618099
E: info@mmsmedical.ie

MMS Medical,
Unit 81,
Western Parkway Business Park,
Lower Ballymount Road,
Dublin 2

T: 01 4500582
E: sales@mmsmedical.ie

www.mmsmedical.ie



BEKA Hospitec GmbH
Am Rübenmorgen 3
35582 Wetzlar-Dutenhofen
info@beka-hospitec.de

Fon +49 641- 92220-0
Fax +49 641- 92220-20
www.beka-hospitec.de

