

Bathroom and Toileting

Everyone wants to retain their own personal freedom – our innovative developments ensure that your customers receive the critical aids they need for their daily care in the bathroom and toilet.

Max. 140 kg



Shown without covers

Made in Germany

KANJO Bathlift



The KANJO bathlift provides safety and freedom during everyday body care. The ergonomically shaped seat and back with adjustable lying positions at the touch of a button guarantee relaxed bathing comfort.

Product details

- For all normal bath tub sizes
- Manual control unit (small) with helical cable and battery charger
- Skin-friendly and non-slip surface
- Lowerable back up to 40°
- Continuously variable raising and lowering at the touch of a button
- Extremely stable base plate (53 x 22 cm)
- Can be folded into 2 main components

Battery

- Lithium-ion technology (charging time: 3 hours)

Manual control with fast charging technology (4 lithium-ion cells)



KANJO

Product description	Bathlift
Seat width without / with side flap (cm)	37.5 / 71
Seat depth (cm)	50
Seat height (cm)	6.5–45.5
Back rest W x H (cm)	35 x 68
Overall height (cm)	110.5
Overall depth, back vertical (cm)	69
Overall depth, back lowered (cm)	90
Max. user weight (kg)	140
Total weight incl. hand control unit (kg)	10.8
Item No.	750/16