



Tilt  
-in-  
Space



Elevating  
Legrest



Back Angle  
Recline



Manual

# FLORIEN ELITE

## Features and Benefits

- Manual tilt-in-space with independent back angle recline
- Manual elevating legrest with negative retraction to aid standing transfers and accommodate contractures
- Removable, "swivel-back" armrests for side transfer and easier hoist access
- Multi-positioning integral headrest
- Adjustable seat width, seat length
- Adjustable seat height (specify at time of manufacture)
- Choice of a sculpted contoured or immersive waterfall backrests
- Ergonomic push handle
- *Ergonomic upholstered arm as standard*
- Simple replacement of modular components
- Reflexion™ foam upholstered in multi-stretch, breathable Dartex on all user contact surfaces
- No Velcro used

## Pressure Relieving Seat Options (integral)

- Flat or ramped removable seat modules:
  - Intelli-Gel® Integral
  - IntelliGent Air
  - Reflexion™ Foam

## Accessories

- Lateral Supports
- Height adjustable, "flip up" Footrest (90° foot positioning) or height adjustable articulating Footrest (Plantar flexion achievable)
- Cushioned or Universal Headrest
- Lap-Belt with D-Ring Straps
- Pelvic Positioners with D-Ring Straps
- Universal Tray

## Weight Limit

- 160kg (353 lbs) / 25 stone

