

Jazzy® Zero Turn

New revolutionary turning technology



iTurn Technology® allows for an incredible **38" turning radius**



Standard features:

- Built in USB charging point
- Front & back monoshock suspension
- LED front & back lights
- Can travel up to 31 miles
- Underseat storage
- Available in 9 different colours
- Front basket

Colour Options



Opal



White



Blue



Black



Pewter



Green



Weight Capacity	Max Speed	Turning Radius	Overall Length	Front/Rear Tyres / Solid tyres	Overall Width	Total Weight w/o batteries	Batteries Required
21 Stone 6 lbs (136kg)	4 Mph (6.4 km/h)	97.2 cm (38.25")	116.2 cm (45.75")	Front: 7" Rear: 9"	69 cm (27")	8 Stone 5 lbs (54 kg)	x2 22Ah Deep Cycle