

# Plinth 2000

## Model 502 2 Section Couch

### Overview:

Robust two-section couch ideally suited for continuous use in areas such as outpatient departments, GP surgeries, colleges, private clinics etc. Featuring smooth anti-entrapment gas operated back rest, lifetime frame and 5 year parts guarantee, this couch forms the backbone of many clinics. A wide range of options are available to expand the use of the couch.



### Standard Features

- UK Manufactured
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Hydraulic or electric operation
- Comfort foam
- Easy operation
- Low running costs
- Hardwearing anti-microbial vinyl
- Safe working load 250kg
- Retractable castors
- Hoist clearance 130mm
- Anti-entrapment back rest

### Overall Dimensions

Length: 188cm (74") approx  
Width: 63cm (25") approx  
Height: 47cm - 97cm approx  
(18" - 38") approx

### Upholstery Dimensions



### Recommended Options:




Motor Back




Battery Back-up



Breathing Hole



Arm Tables




Folding Cotside



Paper Roll Holder

## Operational Parameters



Model 502E - Electric version

Model 502H - Hydraulic version

## For Your Peace of Mind

- Lifetime frame & 5 year parts guarantee (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average lifespan of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2008 by LRQA

## Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl Flame Retardant To CRIB 5 (BS7176 Including Foam)
- Foam To BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

## Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note: Unless specified, beige will be supplied.



MMS Medical,  
51 Eastgate Drive,  
Little Island,  
Cork

T; 021 4618000  
F: 021 4618099  
E : info@mmsmedical.ie  
www.mmsmedical.ie

**Plinth** | 2000

www.plinth2000.com