


Full Back Spacer

The Oxford® Full Back Spacer is an easy to fit contoured sling with padded leg sections, offering good protection for patients who have involuntary movements or behavioural problems. Made from a cushioned and breathable 3D woven 'Spacer' material with high air permeability, offering high levels of comfort.

- ✓ Integrated head support for added support
- ✓ Colour coded binding indicates sling size
- ✓ Positioning handles allow easy patient repositioning
- ✓ Coloured attachment loops promote correct/safe positioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs



Size	Part No.	SWL
XS	SPECIAL	227 kg
S	SL1421	227 kg
M	SL1422	227 kg
L	SL1423	227 kg
XL	SPECIAL	227 kg



Joerns Healthcare strongly recommends that a thorough risk assessment is conducted by a competent person prior to lifting/transferring any client to ensure their complete safety, comfort and dignity is maintained at all times.

MMS Medical, 51 East Gate Drive, Little Island, Cork. Tel: 021 4618000 Email: info@mmsmedical.ie

MMS Medical, Unit 81, Western Parkway business Park, Lower Ballymount Road, Dublin 12. Tel: 01 4500582 Email: sales@mmsmedical.ie

www.mmsmedical.ie