



**MMS Medical**  
51 East Gate Drive,  
Little Island,  
Cork.

**T: 021 461 8000**

**F: 021 461 8099**

**E: [sales@mmsmedical.ie](mailto:sales@mmsmedical.ie)**

**MMS Medical**

Unit 81,  
Western Parkway Business Park,  
Lower Ballymount Road,  
Dublin 12.

**T: 01 450 0582**

**F: 01 456 5889**

**E: [sales@mmsmedical.ie](mailto:sales@mmsmedical.ie)**

# **Easy Reference Product Guide: Cushioning**

## KEY

**Type:** F=Flat P=Postional O=Omni or Bi Directional

**LLL** = Linear Load Limit

**Sizes:** S=Standard C=Custom O=Oversize

**Standard sizes are** 15"x16", 16"x16", 17"x17", 18"x16" and 18"x17"

Standard oversize is 20"x18"

**Weight guide = Recommended weight for a 17"x17" cushion**

*Addition of a 1" crescent base increases the weight guide by approximately 30%.*

*For more weight information please refer to our website – [www.qbitus.co.uk](http://www.qbitus.co.uk).*

Code	Product	Risk Factor	Type	Weight Guide	LLL	Sizes
------	---------	-------------	------	--------------	-----	-------

### COMMUNITY RANGE

GP	General Issue Seat Pad	Low	F & O	14-15st (89-95kg)	21st (133kg)	S, C & O
C1	Community One	Low/Med	F & O	14-15st (89-95kg)	21st (133kg)	S, C & O
C2	Community Two	Med/High	F & O	18st (114kg)	31st (196kg)	S, C & O
C3	Community Three	High	F & O	20st (127kg)	34st (216kg)	S, C & O
CG	Community Gel	Med/High	F & O	14-15st (89-95kg)	21st (133kg)	S & C

### MERCURY RANGE

P.C	Mercury Junior	Low/Med	P	10st (64kg)	21st (133kg)	C
QM50	Mercury 50	Low/Med	P	15st (95kg)	21st (133kg)	S, C & O
QM100	Mercury 100	Low/Med	P	15-16st (95-102kg)	26st (164kg)	S, C & O
QM150	Mercury 150	Med	F & P	17-18st (108-114kg)	31st (196kg)	S, C & O
QM175	Mercury 175	Med	P	17-18st (108-114kg)	31st (196kg)	S, C & O
QM200	Mercury 200	Med/High	P	17-18st (108-114kg)	31st (196kg)	S, C & O
QM300	Mercury 300	High	P	17-18st (108-114kg)	31st (196kg)	S, C & O
QMLP	Mercury Lumbar Pad	—	—	—	—	One Size
QPB	Mercury Posture Back	—	—	—	—	One Size

Code	Product	Risk Factor	Type	Weight Guide	LLL	Sizes
<b>MERCURY XTREME RANGE</b>						
PCX	Mercury Junior Xtreme	Low/Med	P	10st (64kg)	21st (133kg)	C
QM50X	Mercury 50 Xtreme	Low/Med	P	15st (95kg)	21st (133kg)	S, C & O
QM100X	Mercury 100 Xtreme	Low/Med	P	17-18st (108-114kg)	31st (196kg)	S, C & O
QM150X	Mercury 150 Xtreme	Med	F & P	15-16st (95-102kg)	26st (164kg)	S, C & O
QM175X	Mercury 175 Xtreme	Med	P	17-18st (108-114kg)	31st (196kg)	S, C & O
QM200X	Mercury 200 Xtreme	Med/High	P	17-18st (108-114kg)	31st (196kg)	S, C & O
QM300X	Mercury 300 Xtreme	High	P	17-18st (108-114kg)	31st (196kg)	S, C & O
<b>Q-CARE RANGE</b>						
QC1	Q-Care 1	Low/Med	P	16st (102kg)	29st (184kg)	S, C & O
S30	Q-Care 2 (Original)	Med/High	P	16st (102kg)	29st (184kg)	S, C & O
QCINT1	Q-Care Integra 1	Low/Med	P	16st (102kg)	29st (184kg)	S, C & O
QCINT2	Q-Care Integra 2	Med/High	P	18st (114kg)	31st (197kg)	S, C & O
<b>DYNAMIC RANGE</b>						
S9A	Sun-Mate 1"	Low/Med	F	5st (32kg)	—	S, C & O
S9B	Sun-Mate 2"	Med	F	15st (95kg)	29st (184kg)	S, C & O
S9C	Sun-Mate 3"	Med	F	18st (114kg)	31st (197kg)	S, C & O
S10A	Pudgee 1"	Med	F	—	—	S, C & O
S10B	Pudgee 2"	High	F	12st (76kg)	21st (133kg)	S, C & O
S10C	Pudgee 3"	High	F	15st (95kg)	29st (184kg)	S, C & O
S1C	Laminar 2"	High	F	20st (127kg)	36st (229kg)	S, C & O
S1B	Contour	High	F	20-30st (127-191kg)	36st (229kg)	S, C & O
S1A	Super Contour	High	F	20-30st (127-191kg)	36st (229kg)	S, C & O
S78	Interum Cushion	High	P	18st (114kg)	31st (197kg)	S, C & O
S79	Optimus Cushion	High	P	20st (127kg)	31st (197kg)	S, C & O

Code	Product	Risk Factor	Type	Weight Guide	LLL	Sizes
<b>BESPOKE &amp; MISCELLANEOUS</b>						
S43	Castle Top Plus (firm)	Low/Med	F	20st (127kg)	36st (229kg)	S,C & O
S44	Castle Top Plus (medium)	Low/Med	F	16st (102kg)	31st (197kg)	S,C & O
S45	Castle Top (firm)	Low/Med	F	20st (127kg)	31st (197kg)	S,C & O
S46	Castle Top (medium)	Low/Med	F	16st (102kg)	31st (197kg)	S,C & O
S6	Cold Mould 2"	Low/Med	F	13st (83kg)	21st (133kg)	S,C & O
S7	Cold Mould 3"	Low/Med	F	15st (95kg)	23st (146kg)	S,C & O
S8	Cold Mould 4"	Low/Med	F	18st (114kg)	24st (152kg)	S,C & O
S3	Pin-Core 2"	Med	F	8st (51kg)	12st (76kg)	S,C & O
S4	Pin-Core 3"	Med	F	12st (76kg)	14st (89kg)	S,C & O
S5	Pin-Core 4"	Med	F	14st (89kg)	16st (102kg)	S,C & O
S11	Qbi-Flex Cushion 3"	Med/High	F	15st (95kg)	20st (127kg)	S,C & O
S77	Qbi-Gel Cushion	High	F & O	30st (191kg)	34st (216kg)	S & O
DB	Debbonair	Med/High	F & O	20st (127kg)	—	S



# Qbitus Products

MMS Medical,  
51 East Gate Drive,  
Little Island,  
Cork  
T: 021 4618000  
E: [info@mmsmedical.ie](mailto:info@mmsmedical.ie)



[www.mmsmedical.ie](http://www.mmsmedical.ie)