

Tracked Evacuation Chairs

Evacusafe Tracked Evacuation Chairs are designed to enable the safe and comfortable evacuation of mobility impaired people in the event of an emergency, such as fire. The chairs are designed for ease-of-use on internal stairways and fire escapes as well as along flat surfaces.

Our Tracked Evacuation Chairs are suitable for a single person operation. Mainly operated on a straight continuous run, they can also be used (with some training) on a mild spiral, though this will be best suited for a 2 person operation.

Lightweight yet incredibly strong, the chairs can transport persons of up to 182Kg. The tracks create friction which acts as a brake to slow the descent and allows for smooth, controlled operation. Once ground level has been reached the handle can fold away and the chair can then be used as a temporary wheelchair. With a supportive FLAT seat, armrests and padded headrest, the chair is very comfortable to stay in for extended periods.

The chairs fold flat and can be stored on the wall with the hooks supplied along with an 'Evacuation Chair' wall sign. In addition, a tough, long lasting cover is also provided as standard and helps keep the equipment clean and protected ready for use.

We offer two models; the **Excel Evacuation Chair** and the **Standard Evacuation Chair**.



Excel Evacuation Chair

Standard Evacuation Chair

Chair Features

- Ideal for downward stairway evacuation
- Rubberised friction activated tracks
- 4 wheels for optimum manoeuvrability
- Simple set up makes it quick to deploy
- Flat comfortable seat
- Cushioned headrest
- Foot operated brake
- Arm and foot rests
- Built in carry handles
- Quick release safety harness
- Removable covers

	STANDARD CHAIR	EXCEL CHAIR
Weight	12.5Kg	13.6Kg
Load Capacity	150Kg	182Kg
Height (folded/open)	1120mm / 1460mm	1130mm / 1320mm
Depth (folded/open)	150mm / 700 (max)	290mm / 990mm (max)
Width	540mm	650mm
Length in use	1460mm	1430mm