



Raz-AP600

Attendant Propel Heavy Duty

Rehab Shower Commode 600-lb Capacity

The Raz-AP600 [Attendant-Propel, 600-lb (275 kg) capacity] rehab shower commode chair features the Multi-Position Back with fore/aft center of mass adjustment to accommodate a variety of body contours.



MMS
MEDICAL

www.mmsmedical.ie

FEATURES

Extensive range of options available



SPECIFICATIONS	Raz-AP600 (Attendant Propelled 600 lbs capacity)
WEIGHT CAPACITY	600 lbs (275 kg)
SEAT WIDTH (seat and back frame)	22" (56 cm) - 30" (76 cm)
OVERALL WIDTH (with standard armpad)	26.5" (67 cm) with Standard Arm Pad
SEAT WIDTH (commode seat)	23" (58 cm)
SEAT DEPTH (front of back post to front of seat)	18" (46 cm)
OVERALL LENGTH With Z709 Footrests	36" - 39" (91 cm - 99 cm)
CLEARANCE HEIGHT UNDER SEAT	5" (125 mm) casters w/ pan brackets: 14.25" - 18.25" (36 - 46 cm) 5" (125 mm) casters w/o pan brackets: 15" - 19" (38 - 48 cm) 6" (150 mm) casters w/ pan brackets: 15.25" - 19.25" (39 - 49 cm) 6" (150 mm) casters without pan brackets: 16" - 20" (41 - 51 cm)
ARMREST HEIGHT ABOVE SEAT	9.5" (24 cm)
SEAT-TO-FLOOR HEIGHT	19" - 23" in 1" increments (48 - 58 cm)
PRODUCT WEIGHT	5" (125 mm) casters: 54.5 lbs (25 kg) 6" (150 mm) casters: 56.5 lbs (26.7 kg)
FRAME AND HARDWARE MATERIAL	Stainless Steel (Lifetime warranty against corrosion)

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty all other components, subject to normal wear and tear exclusions.

Specifications, and options are subject to change without prior notice.

MANUFACTURED BY RAZ Design Inc. 22 Howden Road, Toronto, Ontario Canada, M1R 3C7

MMS Medical, 51 Eastgate Drive, Little Island, Cork. Tel: 021 4618000 Email: info@mmsmedical.ie

MMS Medical, Unit 81 Western Business Parkway, Lower Ballymount Road, Dublin 12. Tel: 01 4500582

Rev. Date: 2017-02-27