



Raz-SP600

Self Propel Heavy Duty

Rehab Shower Commode 600-lb Capacity

The Raz-SP600 [Self-Propel, 600-lb (275 kg) capacity] rehab shower commode chair features the Multi-Position Back with fore/aft center-of-mass adjustment to accommodate a variety of body contours. Backrest and rear axle plate can be independently adjusted fore/aft to fine-tune center-of-mass position for optimized maneuverability.



MMS
MEDICAL

www.mmsmedical.ie

FEATURES

Extensive range of options available



Tension-adjustable back upholstery

23"W (58 cm) x 18"D (46 cm) contoured, molded seat

Exclusive to Raz, Ischial Pelvic Alignment System (IPAS), provides 2" (5 cm) of fore/aft or 1" (2.5 cm) of lateral seat adjustment

22"W (56 cm) stainless steel frame with 4" (10 cm) of height adjustment

Adjustable, heavy-duty footrests

ALSO AVAILABLE

- Frame widths up to 36" (91 cm)
- Custom frame depths
- Custom seats

SPECIFICATIONS	Raz-SP600 (Self Propelled 600 lbs capacity)
WEIGHT CAPACITY	600 lbs (275 kg)
FRAME WIDTH (seat and back frame)	22" (56 cm) [Custom widths available]
OVERALL WIDTH (with standard armpad)	30" (76 cm) with Standard Arm Pad
SEAT WIDTH (commode seat)	23" (58 cm)
SEAT DEPTH (front of back post to front of seat)	18" (46 cm)
OVERALL LENGTH WHEN UPRIGHT With Z709 Footrests	40" – 42" (101.5 – 106.5 cm)
CLEARANCE HEIGHT UNDER SEAT	5" (125 mm) casters with pan brackets: 15.25" or 16.25" (39 cm or 41 cm) 5" (125 mm) casters without pan brackets: 16" or 17" (40.5 cm or 43 cm) 6" (125 mm) casters with pan brackets: 16.25", 17.25", 20.25" (41 cm, 44 cm, 51.5 cm) 6" (125 mm) casters without pan brackets: 17", 18", or 21" (43 cm, 46 cm, 53 cm)
ARMREST HEIGHT ABOVE SEAT	9.5" (24 cm)
SEAT-TO-FLOOR HEIGHT (5" casters)	20" and 21" (51 cm and 53 cm)
SEAT-TO-FLOOR HEIGHT (6" casters)	20", 21", 24" (51 cm, 53 cm and 61 cm)
PRODUCT WEIGHT	5" (125 mm) casters: 68 lbs (30.9 kg) 6" (150 mm) casters: 69 lbs (31.3 kg)
FRAME AND HARDWARE MATERIAL	Stainless Steel (Lifetime warranty against corrosion)

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions.

Specifications, and options are subject to change without prior notice.

MANUFACTURED BY: RAZ Design Inc. 22 Howden Road, Toronto, Ontario Canada, M1R 3C7

MMS Medical, 51 Eastgate Drive, Little Island, Cork. Tel: 021 4618000 Email: info@mmsmedical.ie

MMS Medical, Unit 81 Western Parkway Business Park, Lower Ballymount Road, Dublin 12. Tel: 01 4500582

Rev. Date: 2017-02-27