



Raz-SP

Self Propel

Rehab Shower Commode

The **RAZ-SP** (Self-Propel) rehab shower commode chair comes standard with 4" (10 cm) of height adjustment, tension-adjustable back upholstery, self-propelling rear wheels, 5" (125 mm) dual-locking casters, flip-up padded armrests, and commode pan. The easy-to-clean, stainless steel frame is designed with unobstructed side and front access for the user or attendant.



MMS
MEDICAL

www.mmsmedical.ie

FEATURES

Extensive range of options available

Tension-adjustable back upholstery

Adjustable back angle

Center-of-mass adjustment

22" (56 cm) or 24" (61 cm) rear wheels

Dual-locking (swivel and roll) casters



Flip-up padded arms

Compatible with all Raz standard-size seats

Exclusive to Raz, Ischial Pelvic Alignment System (IPAS), provides 2" (5 cm) of fore/aft or 1" (2.5 cm) of lateral seat adjustment

Stainless steel frame with 4" (10 cm) of height adjustment

ALSO AVAILABLE

- Long frame – 3" (8 cm) longer for tall clients
- Low frame – 4" seat-to-floor height
- Recliner frame
- 16" and 20" (40 and 51 cm) widths available

SPECIFICATIONS	Raz-SP 22" wheels	Raz-SP 24" wheels
WEIGHT CAPACITY	350 lbs (160 kg)	350 lbs (160 kg)
FRAME WIDTH (seat and back frame)	Standard: 18" (46 cm)	
OVERALL WIDTH (with standard armpad)	26.5" (67 cm)	26.5" (67 cm)
OVERALL LENGTH (with footrest)	39" (99 cm)	39" (99 cm)
SEAT WIDTH (commode seat)	18" (46 cm)	18" (46 cm)
SEAT DEPTH (front of back post to front of seat)	16" or 19" (40.5 or 48 cm)	16" or 19" (40.5 or 48 cm)
CLEARANCE HEIGHT UNDER SEAT	16", 17", 18", 20" (41 cm, 43 cm, 46 cm, 51 cm) without pan 15.25", 16.25", 17.25", 19.25" (39 cm, 41 cm, 44 cm, 49 cm) with pan	15.75", 16.75", 17.75", 19.75" (40 cm, 43 cm, 45 cm, 50 cm) without pan brackets 15", 16", 17", 19" (38 cm, 41 cm, 43 cm, 48 cm) with pan brackets
ARMREST HEIGHT ABOVE SEAT (17" Adjusta-Back)	7" - 10" in ½" increments (18 - 25 cm)	7" - 10" in ½" increments (18 - 25 cm)
SEAT-TO-FLOOR HEIGHT	20" - 24" in 1" increments (51 - 61 cm)	20" - 24" in 1" increments (51 - 61 cm)
PRODUCT WEIGHT (with footrests)	47 lbs (21.3 kg)	49 lbs (22 kg)
FRAME AND HARDWARE MATERIAL	Stainless Steel (Lifetime warranty against corrosion)	Stainless Steel (Lifetime warranty against corrosion)

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions.

Specifications, and options are subject to change without prior notice.

MANUFACTURED BY: RAZ Design Inc. 22 Howden Road, Toronto, Ontario Canada, M1R 3C7 T: (877) 720-5678 E: info@razdesigninc.com

Distributed by : MMS Medical, 51 East Gate Drive, Little Island, Cork. T: 021 4618000 F: 021 4618099 E: info@mmsmedical.ie

MMS Medical, Unit 81 Western Parkway Business Park, Lower Ballymount Road, Dublin 2 T: 01 4500582 E: sales@mmsmedical.ie

www.mmsmedical.ie

RAZSALGO03r01

Rev. Date: 2017-02-27