

PRODUCT DATASHEET

Vomit Bag (Emesis) w/SAP by Curas™ 1.5L, ml-Graduated

Brand	: CURAS [®]
Category	: Nursing
Range	: Medical Consumables
Product Group	: Vomit Bags
REF	: 03051



PRODUCT DESCRIPTION

Disposable PE vomit bag with PVC-free plastic rim, packed with one sachet of 9g SAP by Curas™ (REF 99905). Also known as emesis bag or sick bag. Large 1.5L capacity, with graduation markings to indicate the approximate volume of content.

INTENDED USE

A vomit bag offers a convenient and simple way to manage vomitus and liquid waste (e.g. emesis) in a discreet and hygienic manner. Pre-loaded with a 9g SAP by Curas™ sachet (REF 99905) to solidify liquids and semi-liquids, preventing splashes and eliminating unpleasant odour.

FEATURES & BENEFITS

- Safe and economical alternative to using traditional emesis basins (e.g. kidney dish).
- Wide opening fits over patient's mouth, eliminating spill-over and reducing clean-up costs.
- SAP by Curas™ crystallizes liquid on contact for safer handling of contaminated waste.
- Each vomit bag contains one sachet of SAP by Curas™ for spill and odour prevention.
- Practical, discreet, easy to use and conveniently manage sickness related to gastrointestinal tract, chemotherapy, motion, food poisoning, etc.
- Hard plastic rim makes it easy to hold and handle, thus reducing risk of contamination.
- Strong and deep bag with 1.5L capacity prevents overflow.
- Printed graduation marks for indication of approximate volume.
- "Twist & lock" feature on the rim secures bag content and contains odour for easy disposal.
- Can be used for other liquid waste collection, limiting contamination through leakage.
- Bag printed with simple, illustrated instructions for use.

PRODUCT SPECIFICATIONS

REF	Capacity	Colour	Packing
03051	1.5 L (~51 oz)	White	Inner Bag

Composition

- Vomit bag made of virgin LDPE.
- Rim made of polypropylene material.
- SAP by Curas™ made of superabsorbent polymer (sodium polyacrylate) powder.
- 100% latex-free and PVC-free.

Country of Origin

- China

Sterilization

- Non-sterile

Shelf Life

- 5 years for SAP
- 5 years for Vomit Bag

Selected Standards & Guidelines

- ISO 9001:2008

Recommended Storage Condition

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep in a cool and dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

- If content is not contagious, dispose of closed bag as regular/domestic waste.
- If content is contagious, dispose of closed bag in accordance with clinical waste protocol, or incinerate.

Product Symbols

Vomit Bag

SAP by Curas™



PACKAGING INFO

REF	EAN-13 Barcode	Inner Bag	Carton	Gross Weight
03051	5051400030512	50 pcs / bag	600 pcs / carton (12 inner bags)	~ 14 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
Inner Bag	PE	--	√	(√)		√	
Outer Carton	Corrugated Carton	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√ = Highly recommended

(√) = Recommended

ADDITIONAL INFO

- HS Code: 3926909090
- Documentation available upon request. Kindly contact: info@curas.com